# **BD BACTEC™ MYCO/F-Sputa Culture Vials** Supplemented Middlebrook 7H9 Broth

English: pages 1 – 5 Español: páginas 18 – 22  
Français : pages 5 – 9 Dansk: side 22 – 26  
Deutsch: Seiten 9 – 13 Português: páginas 26 – 30  
Italiano: pagine 14 – 18 Svenska: sidan 31 – 34

PP101JAA(02)  
2015-05

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instrução. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerl BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

## INTENDED USE

**BACTEC™ MYCO/F-Sputa culture media** (modified Middlebrook 7H9 broth with CO<sub>2</sub>) with the addition of **BACTEC™ Supplement/F** and **BACTEC™ PANTA™/F** antibiotic mixture, when appropriate, is used with the **BACTEC™** Brand 9000MB fluorescent series instrument as a qualitative procedure for the *in vitro* culture and recovery of mycobacteria. Acceptable specimens to be used are digested decontaminated clinical specimens and sterile body fluids other than blood.

## SUMMARY AND EXPLANATION

Since the mid-1980s and spread of the AIDS epidemic, *Mycobacterium tuberculosis* (MTB) and *Mycobacteria* other than tuberculosis (MOTT), especially *Mycobacterium avium* complex (MAC), have become resurgent. From 1985 to 1992, the number of MTB cases reported increased 18%. Tuberculosis still kills an estimated 3 million persons a year worldwide, making it the leading infectious disease cause of death.<sup>1</sup> Between 1981 and 1987, AIDS case surveillances indicated that 5.5% of the patients with AIDS had disseminated nontuberculous mycobacterial infections, e.g., MAC. By 1990, the increased cases of disseminated nontuberculous mycobacterial infections had resulted in a cumulative incidence of 7.6%.<sup>2</sup> In addition to the resurgence of MTB, multidrug-resistant MTB (MDR-TB) has become an increasing concern. Laboratory delays in the growth, identification and reporting of these MDR-TB cases contributed at least in part to the spread of the disease.

The Centers for Disease Control and Prevention (CDC) have recommended that every effort must be made for laboratories to use the most rapid methods available for diagnostic mycobacteria testing. These recommendations include the use of both a liquid and a solid medium for mycobacterial culture.<sup>3</sup>

The **BACTEC 9000MB System** is designed for the rapid detection of mycobacteria in clinical specimens other than blood. The system includes a liquid culture medium (MYCO/F-Sputa Culture Vial), a growth supplement (Supplement/F) and an antibiotic supplement (**BACTEC PANTA/F**). **BACTEC Supplement/F** contains growth enhancers for mycobacteria, and is also used to reconstitute **BACTEC PANTA/F**. **BACTEC PANTA/F** contains antimicrobial agents used to suppress the growth of contaminating or normal flora microorganisms which may survive the decontamination process and is recommended as an addition with all non-sterile specimens. Each vial contains a sensor which can detect decreases in the dissolved oxygen in the medium resulting from microorganism metabolism and growth. The sensor is monitored by the instrument for increasing fluorescence which is proportional to decreasing dissolved oxygen content. A positive reading indicates the presumptive presence of viable mycobacteria in the vial. Detection is limited to microorganisms that will grow in the medium.

## PRINCIPLES OF THE PROCEDURE

Clinical specimens other than blood are collected and processed according to standard procedures for the recovery of mycobacteria. Processed specimens are inoculated with a needle and syringe into a **BACTEC MYCO/F-Sputa culture vial** which has been supplemented with **BACTEC Supplement/F** and/or **BACTEC PANTA/F** depending on the specimen type (refer to Procedure section). Patient information is entered into the system computer, and the vial is assigned a station through the bar code menu on the instrument. The vial is placed into the **BACTEC 9000MB System**, and is continuously incubated at 37 °C, with agitation once every ten minutes for maximum recovery. Respiration by viable aerobic microorganisms present in a MYCO/F-Sputa vial results in a net decrease of oxygen within that vial. An increase in the fluorescence of the vial sensor is caused by the depletion of dissolved oxygen. Each culture vial is monitored by the **BACTEC 9000MB System** every ten minutes for increasing fluorescence. Analysis of the rate of dissolved oxygen is used to determine if the vial is culture positive, i.e., the test sample contains viable organisms. Culture vials which remain negative for a minimum of 42 days (up to 56 days), and which show no visible signs of positivity are removed from the instrument and sterilized prior to discarding.

## REAGENTS

The **BACTEC MYCO/F-Sputa culture vials** contain the following active ingredients prior to processing:

List of Ingredients			
Processed Water . . . . .	40 mL	Ammonium Sulfate . . . . .	0.05% w/v
7H9 Middlebrook Broth Base. . . . .	0.47% w/v	Ferric Ammonium Citrate . . . . .	0.006% w/v
Casein Hydrolysate . . . . .	0.10% w/v	Polysorbate 80 . . . . .	0.0025% w/v
Supplement H . . . . .	0.30% w/v	Hemin . . . . .	0.0005% w/v
Glycerol . . . . .	0.10% w/v		

All **BACTEC** media are dispensed with added CO<sub>2</sub>. Composition may have been adjusted to meet specific performance requirements.

Prior to inoculation of respiratory and non-sterile specimens, each 40 mL vial of **BACTEC MYCO/F-Sputa medium** requires the addition of 2.0 mL of **BACTEC PANTA/F** Antibiotic Supplement solution as reconstituted with **BACTEC Supplement/F**. Prior to inoculation of sterile body fluids other than blood, each 40 mL vial of **BACTEC MYCO/F-Sputa medium** must be supplemented with 2.0 mL of **BACTEC Supplement/F**. Please see the Package Inserts for **BACTEC PANTA/F** (PP-102) and **BACTEC Supplement/F** (PP-103) for additional information.

## WARNINGS

**Precautions:** For *in vitro* diagnostic use.  
This Product Contains Dry Natural Rubber.

**POTENTIALLY INFECTIOUS TEST SPECIMEN. Observe “Universal Precautions”<sup>4,5</sup> and institutional guidelines when handling and disposing of infectious materials.**

Biosafety Level 2 practice, containment equipment, and facilities are recommended for preparing acid-fast stains and for culturing clinical specimens. For activities involving the propagation and manipulation of *Mycobacterium tuberculosis* or *Mycobacterium* species grown in culture, Biosafety Level 3 practice, containment equipment, and facilities are required as recommended by the CDC.<sup>6</sup>

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, a vial may exhibit incomplete crimping of the cap, as indicated when the metal edge at the bottom of the cap is not uniformly rolled under the neck finish of the bottle. In both cases the contents of the vials may leak or spill, especially if the vial is inverted. **DO NOT USE** any vial showing evidence of contamination, damage or incomplete crimping.

To minimize the potential of leakage during inoculation of specimen into culture vials, use a tuberculin syringe with a 25-gauge permanently attached needle. A one-handed inoculation technique and a suitable vial holder should be employed to prevent accidental needle stick injury. **Needles larger than 20-gauge must NOT be used.** For supplement addition, pressure release or subculturing, use of larger needles may permanently damage the vial septum, and may cause leakage.

**Positive culture vials for subculturing or staining, etc.:** Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling must be performed in a biological safety cabinet, and appropriate protective clothing, including gloves and masks, should be worn. See **"PROCEDURE"** section for more information on subculturing.

A contaminated vial could contain positive pressure. If a contaminated vial is sampled, gas and/or contaminated culture medium could escape from the vial, creating an aerosol hazard. Vial contamination may not be readily apparent.

Before discarding, sterilize all inoculated MYCO/F-Sputa vials by autoclaving.

#### LEAKING OR BROKEN VIALS

If an inoculated vial is found to be leaking or is accidentally broken, use the established procedure in your laboratory for dealing with mycobacterial spills. As a minimum, "Universal Precautions" should be employed. An approved form of respiratory protection is recommended. Any vial found to have minimal leakage confined to the septum and cap of the vial may be topically disinfected using a mycobactericidal disinfectant, followed by 70% isopropyl alcohol. If this should occur within the instrument, a potential for an aerosol exposure exists. Turn off the instrument and close the doors immediately. Vacate the affected area. Contact your facility's Safety or Infection Control Officer(s). Determine the necessity of turning off or modifying the settings of the air handling units serving the affected area. Do not return to the area until any potential aerosols have settled or have been removed by appropriate ventilation. Becton Dickinson should be notified by calling 1-800-544-7434 in the U.S.A. Guidelines for proper handling of accidental mycobacterial contamination due to breakage of culture tubes or broth suspensions have been issued by the CDC.<sup>6</sup>

#### STORAGE INSTRUCTIONS

**MYCO/F-Sputa culture vials:** Store at 2° – 25 °C, in a dry location out of direct light.

**DO NOT** use after expiration date.

**BACTEC PANTA/F Antibiotic Supplement:** Store unreconstituted **BACTEC PANTA/F** at 2 – 8 °C. Once reconstituted, it should be used immediately or stored frozen at -20 to -70 °C for up to six months. **DO NOT** thaw and refreeze. Protect from light. Overheating must be avoided.

**DO NOT** use after expiration date.

#### SPECIMEN PROCESSING

Sputum or other respiratory specimens must be digested, decontaminated and concentrated prior to inoculation into MYCO/F-Sputa culture vials. The N-acetyl-L-cysteine-sodium hydroxide method is recommended.<sup>6,7,8</sup> Alternatively, the **BBL™ MycoPrep™** Kit may be used for processing the specimen. Other decontamination methods have not been tested in conjunction with MYCO/F-Sputa medium. After centrifugation (≥ 3,000 x g, 15 min) the sediment should be re-suspended in sterile phosphate buffer, pH 6.8.

Processing of non-respiratory specimens other than blood should be conducted according to the Clinical Microbiology Handbook,<sup>9</sup> CDC Manual, or as defined in your individual laboratory procedure manual.

#### PROCEDURE

**MYCO/F-Sputa medium must be used with a BACTEC 9000MB fluorescent series instrument.**

**Materials Provided:** BACTEC MYCO/F-Sputa culture vials, **BACTEC PANTA/F** Antibiotic Supplement and **BACTEC** Supplement/F Media Supplement.

**Materials Required but Not Provided:** Centrifuge; Biological Safety Cabinet; Autoclave; CO<sub>2</sub> incubator at 37 °C, **Falcon™** brand 50 mL centrifuge tubes; 4% sodium hydroxide; 2.9% sodium citrate solution; N-acetyl-L-cysteine powder; phosphate buffer, pH 6.8 vortex mixer; sterile transfer pipettes; sterile tuberculin syringe with 25-gauge needle; mycobactericidal disinfectant; 70% isopropyl alcohol; mycobacterial agar or egg-based medium; 7H9 broth, sterile saline; Quality Control organisms (*Mycobacterium tuberculosis*, ATCC™ 27294; *Mycobacterium fortuitum*, ATCC 6841; and *Mycobacterium kansasii*, ATCC 12478); microscope and materials for staining slides.

#### INOCULATION OF MYCO/F-SPUTA CULTURE VIALS

1. Remove the flip-off cap from the **BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented septum. **DO NOT USE** if any defect is noted.

2. Label culture vial with specimen identification.

3. Reconstitute the lyophilized **BACTEC PANTA/F** vial with 10 mL of **BACTEC** Supplement/F using a syringe with a **Luer-Lok™** brand tip fitted with a 25-gauge needle. Needles larger than 20-gauge **MUST NOT** be used. Larger needles may permanently damage the vial septum. Mix the contents well. Visually inspect the vial contents to insure dissolution of the **BACTEC PANTA/F** Antibiotic Supplement.

4. For respiratory and non-sterile specimens, MYCO/F-Sputa medium vials to be inoculated must be supplemented with 2.0 mL of the reconstituted **BACTEC PANTA/F** solution. Supplemented vials should be inoculated within two hours of the **BACTEC PANTA/F** solution addition. Before inoculating with **BACTEC PANTA/F** solution, swab the septum with 70% isopropyl alcohol. A tuberculin syringe equipped with a 25-gauge, non-removable needle is recommended. A one-handed inoculation technique and a suitable vial holder should be employed to prevent accidental needle stick injury. For tissue or other particulate specimens, it is recommended to use a larger needle (20 – 22 gauge) and a syringe with a **Luer-Lok** brand or similar type tip.

For sterile body fluids other than blood (e.g. synovial fluid, peritoneal), **BACTEC** MYCO/F-Sputa medium vials must be supplemented with 2.0 mL of **BACTEC** Supplement/F prior to inoculation. Before inoculating with **BACTEC** Supplement/F, swab the septum with 70% isopropyl alcohol. A tuberculin syringe equipped with a 25-gauge, non-removable needle is recommended. A one-handed inoculation technique and a suitable vial holder should be employed to prevent accidental needle stick injury.

5. Aseptically inject 0.5 mL of processed specimen per vial using a tuberculin syringe with a 25-gauge needle. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any side-to-side motions which could permanently damage the septum. Swab the septum with mycobactericidal disinfectant, followed by 70% isopropyl alcohol. Inoculated vials should be placed in the **BACTEC** 9000MB System as soon as possible, but certainly within the same day as the specimen is decontaminated, processed, and inoculated into **BACTEC** MYCO/F-Sputa culture vials.

At this time, conventional media such as Lowenstein-Jensen or 7H10/7H11 should be inoculated.

6. Vials entered into the instrument will be automatically tested for the duration of the testing protocol. Positive vials will be identified by the **BACTEC** 9000MB System (see the **BACTEC** 9000MB User's Manual, MA-0092). The sensor inside the vial may not appear visibly different in positive or negative vials; however, the **BACTEC** 9000MB System can determine a difference in sensor fluorescence.

7. Positive vials should be subcultured and an acid-fast smear prepared.

Processing an instrument-positive vial:

- Remove the vial from the instrument.
- In biological safety cabinet, vent the vial to equilibrate vial pressure with atmosphere.
- Invert vial to mix contents.
- Remove aliquot from vial ( app. 0.1 mL) for stain preparations (AFB and Gram stains).
- Inspect smear and preparations. Report preliminary results only after acid fast stain evaluation.

If **AFB positive**, subculture to solid media and report as: instrument-positive, AFB positive, ID pending.

If microorganisms other than acid-fast bacilli are present, report as: instrument positive, AFB Negative, Contaminated.

If no microorganisms are present on the smears, re-enter the vial into the instrument as an ongoing negative vial and allow to complete test protocol. No reportable result.

Subcultures for identification and drug susceptibility testing may be performed using fluid from the **BACTEC MYCO/F-Sputa** culture vials.

**Subculturing and Vial Re-Entry:** Subculturing must be performed in a biological safety cabinet, and appropriate clothing, including gloves and masks, should be worn. Prior to subculturing, place the vial in an upright position, and place an alcohol wipe over the septum. To release any positive pressure in the vial which could be caused by growth of contaminants, insert a sterile 25-gauge (or smaller) needle equipped with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after any pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any side-to-side motions which could permanently damage the septum. To subculture the vented vial, invert to mix contents well and insert a new syringe, with a 25-gauge needle, to remove culture medium for further processing. **Do not re-cap the needle. Discard needles and syringes in a puncture-resistant biohazard container.**

The user may wish to return a vial with a negative smear result to the instrument for continued monitoring through the testing protocol.

At the end of six weeks incubation, perform a visual check of all instrument-negative vials. If the vial appears visually positive (i.e., turbid, with possible clumps of mycobacteria) it should be subcultured, acid-fast stained, and treated as a presumptive positive, provided the acid-fast stain result is positive. If the vial shows no signs of positivity it should be sterilized prior to discarding.

## QUALITY CONTROL

Quality Control Certificates are provided with each carton of media.

It is recommended that each new shipment or lot of **BACTEC MYCO/F-Sputa** culture media be tested with the ATCC control organisms identified in the chart below as a positive control and an uninoculated vial as a negative control.

Organism	Range of Time-to-detection (days)
<i>Mycobacterium tuberculosis</i> , H37Rv, ATCC 27294	8 to 12
<i>Mycobacterium fortuitum</i> , ATCC 6841	1 to 3
<i>Mycobacterium kansasii</i> , ATCC 12478	6 to 12

To prepare the positive control vial:

1. Grow the organism in 7H9 broth.
2. Prepare a #1 McFarland ( $\approx 10^7$  CFU/mL) suspension of the organism using sterile saline.
3. Dilute the suspension to  $10^4$  CFU/mL using sterile saline.
4. Inoculate 1 mL of this suspension into a **BACTEC MYCO/F-Sputa** culture vial which has been supplemented with 2.0 mL of reconstituted **BACTEC PANTA/F**. **Inoculum of 1 mL is to be used for Quality Control organisms only.**

The positive and negative control vials should be scanned into the instrument and tested. The positive control vial should be detected as instrument positive within the range given in the above chart. The negative control should remain negative. If either of these vials do not give the expected results, do not use the media until you have contacted Technical Services at 1-800-638-8663 (United States only).

For information on quality control for the **BACTEC 9000MB** System, refer to the **BACTEC 9000MB** User's Manual (MA-0092).

## RESULTS

A culture-positive sample is determined by the **BACTEC 9000MB** System and confirmed by an acid-fast smear. A positive result indicates the presumptive presence of viable microorganisms in the vial.

## LIMITATIONS OF THE PROCEDURE

### Contamination

Care must be taken to prevent external contamination of the sample during collection, digestion, decontamination, and inoculation into the **BACTEC** vial. Due to the richness of the medium and the non-selective nature of fluorescent detection of oxygen consumption, breakthrough contamination may occur. Digestion and decontamination procedures must be carefully performed to reduce breakthrough contamination by non-mycobacterial organisms. Test only indicated specimen types. A contaminated vial will give a positive instrument reading, but will not indicate a relevant clinical result. Such a determination must be made by the user, based on such factors as acid-fast smear result, type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc. (Refer to the **BACTEC 9000MB** User's Manual, MA-0092, for vial re-decontamination procedure.)

Decontamination with the N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) method is recommended. Other decontamination methods have not been tested in conjunction with **BACTEC MYCO/F-Sputa** medium. Digestant-decontaminants may have harmful effects on mycobacteria.

### Specimen Type

**BACTEC MYCO/F-Sputa** culture vials with the addition of **BACTEC Supplement/F** and **BACTEC PANTA/F**, when appropriate, are to be used for the culture and recovery of digested decontaminated clinical specimens and sterile body fluids other than blood. The minimum detection level of this test may vary according to the mycobacterial species present.

**DO NOT USE BACTEC MYCO/F-Sputa** medium culture vials for the recovery of mycobacteria from blood.

Pleural fluids can contain red and white blood cells<sup>10</sup> which have a high rate of oxidative metabolism which can result in false positivity.

### General Considerations

Detection of mycobacterial species in clinical specimens is dependent on the number of organisms present in the specimen, specimen collection methods, patient factors such as presence of symptoms, prior treatment, and the method of processing. Adherence to procedural instructions is critical for optimum recovery of mycobacteria. Contamination with saprophytic mycobacteria in tap water or other laboratory reagents and equipment may cause false positive results.

Optimum recovery of isolates will be achieved by adding 0.5 mL of processed and re-suspended specimen inoculum to each vial. The severity of the digestion and decontamination protocol used to process the specimen can also affect recovery and time-to-detection.

Mycobacteria may vary in acid-fastness depending on strain, age of culture and other variables. All positive vials as indicated by the instrument or appearing turbid at end of testing protocol should be subcultured to both selective and non-selective mycobacterial media. Non-mycobacterial species may overgrow mycobacteria present. These culture vials should be re-decontaminated and re-cultured. (Refer to the **BACTEC 9000MB** User's Manual, MA-0092, for vial re-decontamination procedure.)

Colonial morphology and pigmentation characteristics can only be determined on solid media. Further identification tests may be performed to determine mycobacterial speciation from a positive **BACTEC MYCO/F-Sputa** vial.

**BACTEC MYCO/F-Sputa** vials which appear positive may contain one or more species of mycobacteria and/or other non-mycobacterial species. Identification of mycobacteria present requires subculture to a solid medium and additional procedures to identify organisms present. The consistency of microscopic morphology in **BACTEC MYCO/F-Sputa** has not been established.


**BACTEC MYCO/F-Sputa** vials are incubated at 37 °C potentially precluding the recovery of mycobacteria requiring other incubation temperatures (e.g., *M. marinum*, *M. ulcerans*, *M. haemophilum*). Recovery of such organisms requires additional culture methods. Organisms with special growth requirements (e.g., *M. haemophilum*) may not be recovered in **BACTEC MYCO/F-Sputa** when incubated at the appropriate temperature. The following isolates were detected as positive in the **BACTEC 9000MB** System using **BACTEC MYCO/F-Sputa** during both internal studies and/or clinical trials: *Mycobacterium abscessus*, *Mycobacterium avium*, *Mycobacterium bovis*, *Mycobacterium celatum*, *Mycobacterium chelonae*, *Mycobacterium fortuitum*, *Mycobacterium goodii*, *Mycobacterium intracellulare*, *Mycobacterium kansasii*, *Mycobacterium malmoense*, *Mycobacterium simiae*, *Mycobacterium szulgai*, *Mycobacterium terrae*, *Mycobacterium tuberculosis* and *Mycobacterium xenopi*.

The use of **BACTEC PANTA/F** Antibiotic Supplement, although necessary for all non-sterile specimens, may have inhibitory effects on some mycobacteria.

Terminal subcultures were not routinely performed during clinical studies. Therefore, an actual false negative rate, defined as a MYCO/F-Sputa culture vial that remained negative throughout the six-week incubation period was subcultured and grew a mycobacterial organism, cannot be determined at this time.

## EXPECTED RESULTS

Frequency distribution of recovery times for clinical trial respiratory specimens positive in the **BACTEC 9000MB** System is illustrated in the following figure.

## PERFORMANCE CHARACTERISTICS

### Respiratory Specimens

The **BACTEC 9000MB** System was evaluated at four clinical sites which included public health laboratories as well as large acute care hospitals in geographically diverse areas. The site populations included patients infected with HIV, immunocompromised patients, and transplant patients. The **BACTEC 9000MB** System was compared to both the **BACTEC 460TB** radiometric system as well as conventional solid growth media for the recovery and detection of mycobacteria from respiratory specimens. A total of 3135 respiratory specimens were tested during the trials. The total number of pathogenic mycobacteria positive isolates recovered in the study was 348. Of these positives, 292 (84%) were recovered in the **BACTEC 9000MB** System, 260 (75%) were recovered in the **BACTEC 460TB** radiometric system and 178 (51%) were recovered using solid media (Lowenstein-Jensen). The **BACTEC 9000MB** System and solid media combined recovered 89.6% of the total pathogenic isolates. For non-pathogenic MOTT\* (Mycobacteria Other Than Tuberculosis), the total number of positive isolates recovered in the study was 40. Of these positives, 16 (40%) were recovered in the **BACTEC 9000MB** System, 24 (60%) were recovered in the **BACTEC 460TB** radiometric system, and five (13%) were recovered in the solid media. The **BACTEC 9000MB** System and solid media combined recovered 50% of the total non-pathogen MOTT.

(\*Non-pathogenic MOTT included *M. goodii*, *M. abscessus*, *M. terrae*.)

The **BACTEC 9000MB** System failed to recover 1.8% of the pathogenic isolates which were recovered in one or more of the reference systems (**BACTEC 460TB** or conventional solid media). While this percentage represents a potential loss of recovery, it is not indicative of an actual false negative determination (refer to Limitations Section). Use of a second medium, as recommended, will increase the probability of recovery of mycobacterial organisms. The **BACTEC 9000MB** System demonstrated a 1.5% false positive rate (instrument-positive, smear and/or subculture negative). The overall breakthrough contamination rate for the **BACTEC 9000MB** System was 6.5%.

### Summary of BACTEC 9000MB System Respiratory isolate recovery during clinical trials

Species Isolated	Total Isolates	TOTAL 9000MB	9000MB ONLY	TOTAL 460TB	460TB ONLY	TOTAL LJ	LJ ONLY
<i>M. tuberculosis</i>	146	123	14	126	13	93	3
<i>M. avium</i> Complex	162	138	38	114	18	65	4
<i>M. kansasii</i>	10	10	0	9	0	9	0
<i>M. fortuitum</i>	18	11	7	6	3	7	4
<i>M. chelonae</i>	5	4	2	3	1	1	0
<i>M. simiae</i>	1	1	1	0	0	0	0
<i>M. malmoense</i>	1	0	0	1	1	0	0
<i>M. goodii</i>	34	14	10	23	19	2	1
<i>M. abscessus</i>	3	2	2	1	1	0	0
<i>M. terrae</i>	1	0	0	0	0	1	1
<i>M. phlei</i>	1	0	0	0	0	1	1
<i>M. vaccae</i>	1	0	0	0	0	1	1
<i>M. species (other)</i>	5	5	2	1	0	3	0
ALL ISOLATES	388	308	76	284	56	183	15

### Non-Respiratory Specimens

In a separate study at a large teaching hospital, 803 non-respiratory specimens were tested with the **BACTEC MYCO/F-Sputa** culture medium, **BACTEC 12B** culture medium, and conventional medium (Lowenstein-Jensen). The total number of pathogenic mycobacteria positive isolates recovered in the study was 38. Of these positive isolates, the **BACTEC 9000MB** System recovered 29 (76.3%), 30 (78.9%) were recovered in the **BACTEC 460TB** System, and 24 (63.2%) were recovered by conventional medium (Lowenstein-Jensen). The **BACTEC 9000MB** System and solid media combined recovered 89.5% of the total pathogenic isolates. Total positive specimens (pathogenic and non-pathogenic mycobacteria) were distributed among the following sources: gastric (5.1%), sterile body fluids other than blood (17.9%), stool (10.3%), superficial skin/wound drainage (5.1%), tissue (53.8%), and urine (7.7%). The following isolates

were detected as positive in the **BACTEC** 9000MB System using MYCO/F-Sputa during this clinical trial: *M. tuberculosis*, *M. avium* complex, *M. chelonae*, *M. fortuitum*, and *M. bovis*. The overall false positive rate (instrument-positive, smear and/or subculture negative) was 5.0%. Due to the variety of specimens collected and tested, the false positive rate varied significantly from the rate previously reported for respiratory specimens. The breakthrough contamination rate for normally sterile specimens (i.e. tissue and sterile body fluids other than blood) ranged from 4.7% – 18.9%; non-sterile specimens (i.e., gastric, stool, urine, superficial skin/wound drainage) ranged from 8.2% – 73.9%. Overall breakthrough contamination rate was 14.9%.

## REFERENCES

1. Bloom, B.R. and Murray, C.J.L., 1992. Tuberculosis: Commentary on a Reemergent Killer. *Science*, 257:1055-1064.
2. Horsburg Jr., C.R. 1991. *Mycobacterium Avium* Complex Infection in the acquired immunodeficiency syndrome. *New England Journal of Medicine* 324:1332-1338.
3. Tenover, F.C., et al., 1993. The resurgence of Tuberculosis: Is Your Laboratory Ready? *Journal of Clinical Microbiology* 31:767-770
4. Recommendations for preventing transmission of Human Immunodeficiency Virus and Hepatitis B Virus to patients during exposure-prone invasive procedures. *MMWR* 1991, Vol. 40, No. RR-8.
5. Bloodborne pathogens. Code of Federal Regulations, Title 29, Part 1910.1030, Federal Register 1991, 56:64175-64182.
6. Kent, P.T. et al. *Public Health Mycobacteriology; A Guide for the Level III Laboratory*, U.S. Department of Health and Human Services Public Health Service/ Centers for Disease Control, Atlanta, GA 30333, 1985. pp. 16-19.
7. Kubica, G.P. et al. Sputum digestion and decontamination with N-acetyl-L-cysteine-sodium hydroxide for culture of mycobacteria. *Amer. Rev. Resp. Dis.* 1963, 87:775-779.
8. Kubica, G.P. et al. Comments on the use of the new mycolytic agent N-acetyl-L-cysteine as a sputum digestant for the isolation of mycobacteria. *Amer. Rev. Resp. Dis.* 1964. 89:284- 286.
9. Master, Ronald N. (ed) *Mycobacteriology*, 1994, Section 3, in Isenberg, Henry D. (ed): *Clinical Microbiology Procedures Handbook*, Vol. 1. Washington, D.C., American Society for Microbiology.
10. Krieg, A.F., 1979. Cerebrospinal fluid and other body fluids, p. 635-679. In J.B. Henry (ed), *Todd-Sanford-Davidsohn, Clinical Diagnosis and Management by Laboratory methods*, 16th ed., vol. 1, W.B. Saunders Co., Philadelphia.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

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Benex Limited
Becton Dickinson Pty Ltd.