# **BD BACTEC™ PANTA™/F Antimicrobial Supplement**

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## INTENDED USE

**BACTEC™ PANTA™/F** Antimicrobial Supplement is recommended as an addition to **BACTEC** Myco/F-Sputa Medium for the primary isolation of mycobacteria from clinical specimens. **PANTA/F** suppresses the normal flora associated with digested and decontaminated specimens. Principal use is with the **BACTEC** brand 9000MB instrument.

## SUMMARY AND EXPLANATION

Since the **BACTEC** system utilizes a liquid medium, it is important to add an antimicrobial supplement to suppress growth of contaminating microorganisms which may survive the decontamination process. **PANTA/F** is a lyophilized mixture of five antimicrobials that is reconstituted with Supplement/F. Two mL of the reconstituted **PANTA/F** is added to each Myco/F-Sputa Medium vial which then is inoculated with a specimen. The culture vial is incubated and tested on a **BACTEC** brand 9000MB instrument for the detection of mycobacteria.

## PRINCIPLES OF THE PROCEDURE

The **PANTA/F** product contains polymyxin B, amphotericin B, nalidixic acid, trimethoprim and azlocillin, which are antimicrobial agents used to suppress the growth of contaminating or normal flora microorganisms which may survive the decontamination process. **PANTA/F** is reconstituted with Supplement/F and then added to Myco/F-Sputa Medium. **PANTA/F** is used with specimens expected to contain contaminating microorganisms. However, Supplement/F is a required component of the system regardless of whether contamination is suspected [see the Supplement/F package insert (PP103JAA)].

## REAGENTS

**PANTA/F** vials contain the following active ingredients prior to processing:

Ingredients per vial of lyophilized agents:

Approximate formula per vial:	
Polymyxin B . . . . .	10,000 units
Amphotericin B . . . . .	1,000 µg
Nalidixic Acid . . . . .	4,000 µg
Trimethoprim . . . . .	1,000 µg
Azlocillin . . . . .	1,000 µg

## Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. “Standard Precautions”<sup>1-4</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.**

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials that exhibit evidence of damage should be discarded prior to inoculation. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill.

For use with Myco/F-Sputa Medium; read and follow instructions and information (PP101JAA) supplied with that product.

Reconstitute the **PANTA/F** vial only with Supplement/F. DO NOT add **PANTA/F** to Myco/F-Sputa Medium culture vials *more than 2 h* prior to the inoculation of specimens.

## Storage Instructions

Store unreconstituted **PANTA/F** at 2 – 8 °C. Once reconstituted, it should be used immediately or stored frozen at -20° to -70 °C for up to 6 months. **Do not** thaw and refreeze. Protect from light. Overheating must be avoided.

Do not use after expiration date.

## SPECIMEN COLLECTION AND PREPARATION

See the Myco/F-Sputa Medium package insert (PP101JAA) for directions.

## PROCEDURE

**Materials Provided:** BACTEC PANTA/F Antimicrobial Supplement

**Materials Required But Not Provided:** BACTEC Myco/F-Sputa Medium, BACTEC Supplement/F Media Supplement, BACTEC 9000MB Instrument.

CAUTION: Examine vials for evidence of contamination or defects. Do not use any defective or contaminated vials.

### Preparation of PANTA/F Supplement:

Remove the metal cap and rubber septum from a PANTA/F vial. Aseptically transfer 10.0 mL of BACTEC Supplement/F into the vial. Gently mix until completely dissolved.

### Addition of PANTA/F Supplement to BACTEC Myco/F-Sputa Medium:

Before inoculating, swab the septum with mycobactericidal disinfectant, followed by 70% alcohol. Aseptically add 2.0 mL of reconstituted PANTA/F to the Myco/F-Sputa Medium vial not more than 2 h prior to inoculation of the specimen. A tuberculin syringe equipped with a 25-gauge, non-removable needle or Luer-Lok™ brand tip is recommended. A one-handed inoculation technique and a suitable vial holder should be employed to reduce the risk of accidental needle stick injury. Needles larger than 20-gauge should not be used, as larger needles may permanently damage the vial septum.

## QUALITY CONTROL

Do not use vials past their expiration date.

DO NOT use vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each kit of BACTEC PANTA/F. Quality control certificates show test organisms utilized for this type of product.

For information on quality control for the BACTEC 9000MB instrument, refer to the BACTEC 9000MB User's Manual (MA-0092). Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

See the Myco/F-Sputa Medium package insert PP101JAA.

## LIMITATIONS OF THE PROCEDURE

PANTA/F has been designed to minimize overgrowth of contaminating flora; however, when sputum specimens have been insufficiently or improperly decontaminated, overgrowth can occur. Contaminants resistant to the antimicrobials contained in PANTA/F may be encountered.

## EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

See the Myco/F-Sputa Medium package insert PP101JAA.

## AVAILABILITY

Cat. No.	Description
442188	BD BACTEC™ PANTA™/F Antimicrobial Supplement, Pkg. of 10 vials
442189	BD BACTEC™ Supplement/F Media Supplement, Pkg. of 10 vials
442187	BD BACTEC™ Myco/F-Sputa Medium, Ctn. of 50 vials

## REFERENCES

1. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or [www.bd.com/ds](http://www.bd.com/ds).

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Benex Limited
Becton Dickinson Pty Ltd.