BD BACTEC™ Supplement/F Media Supplement

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INTENDED USE

BACTEC™ Supplement/F Media Supplement is a required supplement to Myco/F-Sputa Medium for the recovery of mycobacteria. In addition, it is used to reconstitute **PANTA™/F** Antimicrobial Supplement. Principal use is with Myco/F-Sputa Medium and with the **BACTEC** brand 9000MB instrument.

SUMMARY AND EXPLANATION

Supplement/F is used to reconstitute **PANTA/F** Antimicrobial Supplement. The reconstituted fluid is then added to each Myco/F-Sputa Medium vial prior to specimen inoculation. Supplement/F contains essential nutrients for the growth and detection of mycobacteria and therefore is a required supplement to Myco/F-Sputa Medium even if **PANTA/F** antibiotics are not needed. The Myco/F-Sputa Medium culture vial is incubated and tested in a **BACTEC** 9000MB instrument for the detection of mycobacteria.

PRINCIPLES OF THE PROCEDURE

Supplement/F contains polyoxyethylene stearate (POES), bovine serum albumin and other nutrients which are growth enhancers for mycobacteria. Supplement/F is also used to reconstitute **PANTA/F**.

REAGENTS

Supplement/F vials contain the following active ingredients prior to processing:

List of Ingredients

Purified Water	12 mL
Lactic Acid	3.5 mg/mL
Polyoxyethylene Stearate	2.3 mg/mL
Bovine Serum Albumin	116.0 mg/mL
Dextrose	23.0 mg/mL
Biotin	0.012 mg/mL

Warnings and Precautions

For *in vitro* diagnostic use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials that exhibit evidence of damage should be discarded prior to inoculation. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill.

For use with Myco/F-Sputa Medium; read and follow instructions and information supplied with that product.

Reconstitute the **PANTA/F** vial only with Supplement/F. **DO NOT** add the supplemented fluid to Myco/F-Sputa Medium culture vials *more than 2 h* prior to the inoculation of the specimen.

Storage Instructions

Store at 15 - 25 °C. See **PANTA/F** package insert (PP102JAA) for storage instructions after the Supplement/F has been used to reconstitute **PANTA/F** Supplement.

Do not use after expiration date.

SPECIMEN COLLECTION AND PREPARATION

See the Myco/F-Sputa Medium package insert (PP101JAA) for directions.

PROCEDURE

Materials Provided: **BACTEC** Supplement/F Media Supplement

Materials Required But Not Provided: BACTEC Myco/F-Sputa Medium, BACTEC PANTA/F Antimicrobial Supplement, BACTEC 9000MB Instrument.

CAUTION: Examine vials for evidence of contamination or defects. Do not use any defective or contaminated vials.

Reconstitution of PANTA/F Antimicrobial Supplement:

Remove the metal cap and rubber septum from a PANTA/F vial. Aseptically transfer 10.0 mL of Supplement/F into the vial. Gently mix until completely dissolved.

Addition of PANTA/F to BACTEC™ Myco/F-Sputa Medium:

See the PANTA/F Supplement package insert (PP102JAA) for directions.

Addition to BACTEC Myco/F-Sputa Medium without PANTA/F:

Supplement/F is a required component of the system. If PANTA/F is not required (as, for example, with QC organisms), Supplement/F *must* still be added prior to inoculation of specimen or organism suspension. Before inoculating, swab the septum with 70% alcohol and allow to dry. Aseptically add 2.0 mL of Supplement/F to the Myco/F Medium vial. A tuberculin syringe equipped with a 25-gauge, non-removable needle or Luer-Lok™ brand tip is recommended. A one-handed inoculation technique and a suitable vial holder should be employed to reduce the risk of accidental needle stick injury. Needles larger than 20-gauge should not be used, as larger needles may permanently damage the vial septum.

QUALITY CONTROL

Do not use vials past their expiration date.

DO NOT use vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each kit of BACTEC Supplement/F media. Quality Control Certificates show test organisms utilized for this type of product.

For information on quality control for the BACTEC 9000MB instrument, refer to the BACTEC 9000MB User's Manual (MA-0092).

See BACTEC Myco/F-Sputa Medium package insert, (PP101JAA), for medium quality control.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

See BACTEC Myco/F-Sputa Medium package insert PP101JAA.

LIMITATIONS OF THE PROCEDURE

See BACTEC Myco/F-Sputa Medium package insert PP101JAA.

EXPECTED RESULTS AND PERFORMANCE CHARACTERISTICS

See BACTEC Myco/F-Sputa Medium package insert PP101JAA.

AVAILABILITY

Cat. No. Description

442189	BD BACTEC™ Supplement/F Media Supplement, Pkg. of 10 vials
442188	BD BACTEC™ PANTA™/F Antimicrobial Supplement, Pkg. of 10 vials
442187	BD BACTEC™ Myco/F-Sputa Medium, Ctn. of 50 vials

REFERENCES

1. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
3. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD BACTEC Supplement/F Media Supplement

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Benex Limited
Becton Dickinson Pty Ltd.