
Revisions

Rev from	Rev to	ECO #
0108	2010/02	5262-10

Notes:

- BD Cat. Number 442206
- Blank (Sheet) Size : Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): # N/A

- See Specification Control Number VS-13-7240-0PR for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: N/A PMS# N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP113JAA		Category and Description Package Insert, BACTEC Mycosis-IC/F	Sheet: 1 of 25 Scale: N/A	A

BD BACTEC™ Mycosis-IC/F Culture Vials

Selective Medium for Yeast and Fungi

English:	pages	1 – 3	Español:	páginas	10 – 12
Français :	pages	3 – 5	Dansk:	side	12 – 14
Deutsch:	Seiten	5 – 7	Português:	páginas	14 – 16
Italiano:	pagine	8 – 10	Svenska:	sidan	16 – 18

 PP113JAA
2010/02

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свяжитесь с местним представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндіңіз жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BACTEC™ Mycosis-IC/F culture vials are for aerobic blood cultures. Principal use is with the **BACTEC** fluorescent series instruments for the selective culture and recovery of yeasts and fungi from blood.

SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the **BACTEC** fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO₂ produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO₂ present. A positive reading indicates the presumptive presence of viable microorganisms in the vial.

PRINCIPLES OF THE PROCEDURE

If yeast and fungi are present in the test sample inoculated into the **BACTEC** vial, CO₂ will be produced when the microorganisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the increase of CO₂ are monitored by the **BACTEC** fluorescent series instrument. Analysis of the rate and amount of CO₂ increase enables the **BACTEC** fluorescent series instrument to determine if the vial is positive, i.e., that the test sample contains viable organisms.

REAGENTS

The **BACTEC** Mycosis-IC/F culture vials contain the following reactive ingredients prior to processing:

List of Ingredients

Processed Water	40 mL
Brain Heart Infusion Broth	1.0% w/v
Soybean-Casein Digest Broth	0.5% w/v
Yeast Extract	0.035% w/v
Sucrose	0.6% w/v
Dextrose	0.1% w/v
m-Inositol	0.05% w/v
Ferric Ammonium Citrate	0.0001% w/v
Sodium Polyanetholesulfonate (SPS)	0.05% w/v
Saponin	0.24% w/v
Chloramphenicol	0.0037% w/v
Tobramycin	0.001% w/v
Antifoaming Agent	0.01% w/v

All **BACTEC** media are dispensed with added CO₂. Composition may have been adjusted to meet specific performance requirements.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of damage, contamination or deterioration. Vials displaying evidence of damage or contamination such as leakage, cloudiness, discoloration (darkening), bulging or depressed septum should not be used.

A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. When using direct draw procedures, monitor the process closely to avoid refluxing materials into patient.

On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vial may leak or spill. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or Luer-Lok™ brand tips.

Storage Instructions

The **BACTEC** vials are ready for use as received and require no reconstitution or dilution. Store in a dry place at 2–25°C, out of direct light.

SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The typical specimen volume is 8 – 10 mL. It is recommended that the specimen be inoculated into the BACTEC vials at bedside. Most commonly, a 10 cc or 20 cc syringe with a Luer-Lok brand tip is used to draw the sample. If appropriate, a Vacutainer™ Brand Needle Holder and a Vacutainer™ Brand Blood Collection Set, Vacutainer™ Safety-Lok™ Blood Collection Set or other tubing “butterfly” set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 8 – 10 mL has been drawn, the flow should be stopped by crimping the tubing and removing the tubing set from the BACTEC vial. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. **The inoculated BACTEC vial should be transported as quickly as possible to the laboratory.** A yellow-top Vacutainer™ Brand Tube containing SPS may also be used to collect the blood sample from the patient. The tube should be transported to the laboratory as quickly as possible for transfer into the BACTEC culture vial.

PROCEDURE

Remove the flip-off cap from BACTEC vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented stoppers. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is not recommended). Aseptically inject or draw directly 8 – 10 mL of specimen per vial. The medium is designed for use with blood specimens ranging from 3 to 10 mL. If the sample is less than 3 mL, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated vials should be placed in the BACTEC fluorescent series instrument as soon as possible for incubation and monitoring.** If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the BACTEC fluorescent series instrument, but rather it should be subcultured, appropriately stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the BACTEC fluorescent series instrument and identified as such (see the appropriate BACTEC fluorescent series instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the BACTEC fluorescent series instrument can determine a difference in fluorescence.

If at the end of the testing period a negative BACTEC Mycosis-IC/F vial appears visually positive (i.e., bulging septum, and/or turbid), it should be subcultured, appropriately stained and treated as a presumptive positive.

Positive vials should be subcultured and stained appropriately. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician.

Subculturing: Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

For maximum yield of isolates, negative cultures may be checked by stain and/or subcultured at some point prior to discarding as negative.

QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

DO NOT USE culture vials past their expiration date.

DO NOT USE culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates show test organisms, including ATCC™ cultures specified in the CLSI Standard, *Quality Control for Commercially Prepared Microbiological Culture Media*.⁵

The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

Organism

*Candida albicans** ATCC 10231

Candida parapsilosis ATCC 10232

Candida (Torulopsis) glabrata ATCC 15545

Candida albicans ATCC 14053

Candida krusei ATCC 34135

Cryptococcus neoformans ATCC 13690

Candida tropicalis ATCC 750

*CLSI Strain

It is recommended that each shipment of media be tested for performance through the use of a positive and negative vial test. The positive vial should be inoculated with 0.1 mL of a 0.5 McFarland Standard of either *Candida albicans* (ATCC 10231) or *Candida (Torulopsis) glabrata*. This vial and an uninoculated vial should be logged into the instrument and tested. The inoculated vial should be detected as positive by the instrument within 72 hours. The negative control vial should remain negative. This verifies that the media were not subject to adverse storage or shipping conditions prior to receipt in your laboratory. If either of these vials do not give the expected results, do not use the media until you have contacted Technical Services at 1-800-638-8663, or your BD representative.

For information on Quality Control for the BACTEC fluorescent series instrument, refer to the appropriate BACTEC fluorescent series instrument User's Manual.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the BACTEC vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organisms recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Nonviable Organisms

A Gram-stained smear from culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.⁶

Dimorphic fungi

Mycosis-IC/F Medium has been optimized for the recovery of yeast. The medium has not been proven effective for the recovery of dimorphic fungi.

General Considerations

Optimum recovery of isolates will be achieved by adding 8 – 10 mL of blood. Use of lower or higher volumes may adversely affect recovery and/or detection times. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO₂ to be detected by the system or significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high.

Due to the nature of biological materials in media products and inherent organism variability, the user should be cognizant of potential variable results in the recovery of certain microorganisms.

EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

Performance in Seeded Blood Culture Studies

Seeded blood culture studies were performed using whole human blood and inoculum levels of 10 to 100 CFU. The following is a list of the organisms which grew in Mycosis-IC/F medium and were detected in the BACTEC fluorescent series instruments.

List of Yeast & Fungi Detected in BACTEC Mycosis-IC/F Medium in Seeded Blood Culture Studies

<i>Aspergillus flavus</i>	<i>Candida parapsilosis</i>	<i>Mucor ramosissimus</i>
<i>Candida albicans</i>	<i>Candida rugosa</i>	<i>Penicillium roquefortii</i>
<i>Candida glabrata</i>	<i>Candida tropicalis</i>	<i>Trichosporon capitatum</i>
<i>Candida krusei</i>	<i>Cryptococcus neoformans</i>	

The performance of BACTEC medium for the recovery of yeast and fungi has been previously established with BACTEC NR Fungal medium and BACTEC non-radiometric systems. Seeded laboratory studies performed by BD have shown equivalent performance of the Mycosis-IC/F medium to NR Fungal medium.

AVAILABILITY

Cat. No.	Description
442206	BACTEC™ Mycosis-IC/F Culture Vials, case of 50 vials

REFERENCES

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17: 53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control for commercially prepared microbiological culture media, 3rd ed., CLSI, Wayne, Pa.
6. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.). 2007. Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.

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Benex Limited