
**Revisions**

Rev from	Rev to	ECO #
0607	2010/07	5394-10

**Notes:**

- BD Cat. Number 442104
- Blank (Sheet) Size : Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number VS-0250 for Material Information
- Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# 199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP117JAA		Category and Description Package Insert, BACTEC Diluting Fluid	Sheet: 1 of 17 Scale: 1:1	<b>A</b>

# **BD BACTEC™ Diluting Fluid**

English: pages 1 – 3 Italiano: pagine 7 – 8  
Français : pages 3 – 5 Portuguais: páginas 8 – 10  
Deutsch: Seiten 5 – 6 Español: páginas 10 – 12

PP117JAA  
2010/07

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свържете се с местния представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / **Обратитесь к своему локальному представителю компании BD за инструкцией.** / Для получения инструкций свяжитесь с местным представителем компании BD. / **Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз.** / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

**BACTEC™ Diluting Fluid** is used in procedures where mycobacterial cultures are diluted for inoculation. Its main use is in the standard **BACTEC** drug susceptibility testing of *M. tuberculosis*.<sup>1,3-7</sup> In this procedure the inoculum is diluted 1:100 for the control vial. Diluting Fluid may also be used to dilute cultures in other procedures such as the NAP test or any other test where dilution of a test culture suspension is required.<sup>2</sup>

## SUMMARY AND EXPLANATION

Use of Diluting Fluid is a part of the **BACTEC** drug susceptibility testing of *M. tuberculosis*. The Diluting Fluid contains a detergent, polysorbate 80 (Tween 80), and a fatty acid-free albumin. These substances help in the dispersal of mycobacteria in a suspension. The albumin is fatty-acid free so that it does not compete with <sup>14</sup>C labelled fatty acid present in the medium. The concentration of polysorbate 80 is minimal so that it may not interfere in the drug susceptibility test but still helps in the dispersion of mycobacteria to give a uniform and homogenous suspension.

## PRINCIPLES OF THE PROCEDURE

The **BACTEC** Susceptibility Test is based on a modified proportion method. The critical proportion for resistance is taken as 1% which means that if 1% or more of the test mycobacterial population is resistant, the culture is considered resistant. To determine the 1% proportion of resistance in 12B Mycobacteria Medium, the bacterial inoculum for the control vial is first diluted 1:100 before adding it to the medium vial. **BACTEC Diluting Fluid** is provided to facilitate the 1:100 dilution of the culture suspension. Fatty acid-free albumin and polysorbate 80, present in the Diluting Fluid, help in dispersing the suspension and give a more uniform suspension of mycobacteria.

## REAGENTS

The **BACTEC Diluting Fluid** vials contain the following reactive ingredients prior to processing:

### **BACTEC DILUTING FLUID, 9.9 ML EACH**

Amount added per Liter of Processed Water

Polysorbate 80.....	0.2 mL
Fatty Acid-Free Bovine Serum Albumin.....	2.0 g

## Warnings and Precautions:

**BACTEC Diluting Fluid** is for *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>8-11</sup> and institutional guidelines should be followed in handling all items contaminated with blood or other body fluids.

Prior to use, each vial should be examined for evidence of damage, deterioration or contamination such as cloudiness, bulging or depressed septum. **DO NOT USE** any vial showing evidence of contamination. A contaminated vial could contain positive pressure. Vial contamination may not be readily apparent. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. Vials that exhibit evidence of damage should be discarded prior to inoculation. On rare occasions, the glass bottle may be cracked and the neck may break during removal of the flip-off cap or in handling. Also on rare occasions a vial may not be sealed sufficiently. In both cases, the contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Aseptic techniques and established precautions against microbial hazards should be used at all times. Procedures should be performed in a suitable biological safety cabinet in a room with an appropriate ventilation system as recommended by the CDC.<sup>12</sup> Use a proper protective gown, mask, and gloves while handling specimens and cultures of potential pathogens. Follow the CDC and OSHA recommendations.<sup>12</sup>

To minimize the potential of leakage during inoculation, use syringes with permanently attached needles or securely fastened **Luer-Lok™** brand tips.

## STORAGE INSTRUCTIONS

**BACTEC** Diluting Fluid vials should be stored at 2 – 25°C, in a dry place out of direct sunlight.

## SPECIMEN COLLECTION

See instructions for susceptibility testing (S.I.R.E. Drug Kit package insert, PP-118JAA).

## PROCEDURE

### Materials Provided:

**BACTEC** Diluting Fluid

### Materials Required But Not Provided:

**BACTEC** 12B Medium

**BACTEC** 460TB Instrument

Tuberculin syringe with permanently attached needle

Disinfectant solution

70% isopropyl alcohol

Remove the flip-off cap from the **BACTEC** Diluting Fluid vial. The rubber septum of a vial should be disinfected each time the vial is entered. Clean the top of each vial with disinfectant, followed by 70% alcohol. Add the culture suspension to be diluted to the vial using a syringe with a permanently attached needle or securely fastened **Luer-Lok** brand tip. When 0.1 mL is added to the vial, a 1:100 dilution is achieved. Mix thoroughly by inverting a minimum of 10 times. If dilutions other than 1:100 are required, remove an aliquot from the vial and make dilutions as appropriate. For further details of the drug susceptibility test procedure, see the S.I.R.E. Drug Kit package insert (PP-118JAA) or the **BACTEC** 460TB Product and Procedure Manual (MA-0029).

## QUALITY CONTROL

DO NOT USE vials past their expiration date.

DO NOT USE vials that exhibit any cracks or defects; discard the vial in an appropriate manner.

The user should examine the Diluting Fluid vials for evidence of deterioration prior to use. For a quality control check of bacterial procedures, please refer to the S.I.R.E. Drug Kit package insert (PP-118JAA) or the **BACTEC** 460TB System Product and Procedure Manual (MA-0029).

Quality Control Certificates are provided with each Diluting Fluid Kit. Quality Control Certificates show test organisms utilized for this type of product.

As part of the daily maintenance schedule, Performance Test Vials (**BACTEC** Performance Test Kit) should be tested on the **BACTEC** 460TB System instrument. If low readings are obtained, yet Performance Test Vials were tested properly, possible instrument problems may be indicated (see the **BACTEC** PTK package insert, PP-046JAA).

## LIMITATIONS OF THE PROCEDURE

The accurate preparation of a 1:100 dilution is achieved only when the suspension to be diluted has been thoroughly dispersed. The presence of large clumps, especially in the case of *M. tuberculosis*, may result in incorrect dilution and thus may affect the drug susceptibility results.

### Contamination

Care must be taken to prevent contamination during inoculation into the **BACTEC** vial. A contaminated sample will give a positive reading but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, appropriate equipment checks, etc. Contamination may also result from improper maintenance of **BACTEC** needles and tubing, poor septum cleansing or malfunction of the **BACTEC** needle heating unit (see Operation and Maintenance Manual).

Due to the nature of biological materials in media products and inherent organism variability, the user shall be cognizant of potential variable results in the recovery of certain organisms.

## EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

See the S.I.R.E. Drug Kit package insert (PP-118JAA) for details.

When a culture suspension such as *M. tuberculosis* is mixed in the Diluting Fluid after thoroughly mixing a homogenous suspension of the culture suspension, a desired dilution is achieved. The presence of polysorbate 80 does not significantly affect the susceptibility pattern due to its low concentration.

## REFERENCES

1. Hawkins, J.E. Non-weekend schedule for **BACTEC** susceptibility testing of *Mycobacterium tuberculosis*. J. Clin. Microbiol. 23:934-937, 1986.
2. Laszlo, A. and S.H. Siddiqi. Evaluation of a rapid radiometric differentiation test for the *Mycobacterium tuberculosis* complex by selective inhibition with r-nitro-a-acetylamino-b-hydroxy-propiofenone. J. Clin. Microbiol. 19:694-695, 1984.
3. Libonati, J.P. et al. Identification and direct drug susceptibility testing of mycobacteria by the radiometric method. Abstract C300, A.S.M. Annual Meeting, St. Louis, MO, 1984.
4. Roberts, G.D. et al. Evaluation of **BACTEC** radiometric method for recovery of mycobacteria and drug susceptibility testing of *Mycobacterium tuberculosis* from acid-fast smear-positive specimens. J. Clin. Microbiol. 18:689-696, 1983.
5. Siddiqi, S.H. et al. Evaluation of a rapid radiometric method for drug susceptibility testing of *Mycobacterium tuberculosis*. J. Clin. Microbiol. 1981, 13:908-912.

- Siddiqi, S.H. et al. Interlaboratory drug susceptibility testing of *Mycobacterium tuberculosis* by radiometric and two conventional methods. *J. Clin. Microbiol.* 1985, 22:919-923.
- Snider, Jr., D.E. et al. Rapid drug susceptibility testing of *Mycobacterium tuberculosis*. *Am. Rev. Respir. Dis.* 123:402-406, 1981.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed., NCCLS, Wayne, PA.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
- Kent, P.T. et al. Public health mycobacteriology, a guide for the level III laboratory. Centers for Disease Control, Division of Laboratory Training and Consultation. Atlanta, GA, 1985.

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Benex Limited