
**Revisions**

Rev from	Rev to	ECO #
E	F	5359-10

**Notes:**

- BD Cat. Number 240952, 240953
- Blank (Sheet) Size: Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): N/A

- See Specification Control Number 240952 for Material Information
- Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS# N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: X5401		Category and Description Package Insert, BBL Staphyloslide Latex Test	Sheet: 1 of 31 <hr/> Scale: N/A	<b>A</b>

# **BD BBL™ Staphyloslide™ Latex Test** for *Staphylococcus aureus*

English: pages 1 – 6 Italiano: pagine 17 – 21  
Français : pages 6 – 11 Español: páginas 22 – 26  
Deutsch: Seiten 11 – 16

X5401F  
2010/06

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrucție ziskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

## INTENDED USE

The **BBL™ Staphyloslide™ Latex Test** is a latex slide agglutination test for the differentiation of staphylococci which possess clumping factor and/or Protein A, usually present with *Staphylococcus aureus*, from staphylococci that do not possess these properties.<sup>1</sup>

## SUMMARY AND EXPLANATION

It has been reported that approximately 97% of human strains of *S. aureus* possess both bound coagulase and extracellular staphylocoagulase. Protein A is found on the cell surface of about 95% of human strains of *S. aureus* and has the ability to bind the Fc portion of immunoglobulin G (IgG).<sup>2</sup>

## PRINCIPLES OF THE PROCEDURE

Traditionally, differentiation between coagulase positive and negative staphylococci has been performed with the tube coagulase test that detects extracellular staphylocoagulase or the slide coagulase test that detects the clumping factor (bound coagulase) present on the bacterial cell surface. Several other differentiation tests are also available, including the passive hemagglutination test and the DNase test.

The **BBL Staphyloslide** Latex Test consists of blue latex particles coated with human fibrinogen and IgG. On mixing the latex reagent with colonies of staphylococci which have clumping factor or Protein A present, cross linking will occur giving visible agglutination of the latex particles. Such agglutination will occur notably with *S. aureus*. If neither clumping factor nor Protein A are present, no agglutination will occur and the result will be regarded as negative. The most

frequent coagulase and Protein A negative isolates of staphylococci are *Staphylococcus epidermidis*.

REAGENTS		100 TESTS	500 TESTS
<b>BBL Staphyloslide Test Latex</b>	Blue latex particles coated with human fibrinogen and IgG, with 0.1% sodium azide (preservative).	2 x 2.5 mL	10 x 2.5 mL
<b>BBL Staphyloslide Control Latex</b>	Blue unsensitized latex particles, with 0.1% sodium azide (preservative).	2 x 2.5 mL	10 x 2.5 mL
<b>BBL Staphyloslide Latex Reaction Cards</b>	Disposable; each card may be used for 6 tests.	35 cards	175 cards

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe established precautions against microbiological hazards throughout all procedures. After use, contaminated materials must be sterilized by autoclaving.

**REAGENTS:** The test latex can become contaminated if the dropper tip is allowed to touch the specimen on the reaction card. Ensure caps on reagent bottles are securely fitted after each use to prevent contamination and drying out of the reagents.

**WARNING:** The reagents contain material of human origin. The human plasma proteins used in the manufacture of the reagent have been tested for the presence of the antibody to HIV (Human Immunodeficiency Virus), Hepatitis C and HBsAg (Hepatitis B surface Antigen) and found to be nonreactive. Because no test method can offer complete assurance that infectious agents are absent, **SPECIMENS AND THESE REAGENTS SHOULD BE HANDLED AS THOUGH CAPABLE OF TRANSMITTING AN INFECTIOUS DISEASE.** The U.S. Food and Drug Administration (FDA) recommends such material be handled at a Biosafety Level 2. BSL 2 is referenced in the U.S. Centers for Disease Control and Prevention/National Institutes of Health (CDC/NIH) manual, *Biosafety in Microbiological and Biomedical Laboratories*.

Reagents contain sodium azide, which is very toxic by inhalation, in contact with skin, and if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**CARDS:** Cards must be flat for proper reactions. If necessary, flatten cards by bowing back in a direction opposite to that of the curl. Care should be taken not to finger-mark the test areas, since this may result in an oily deposit and improper test result. Use each card once and discard. Store cards in the original package in a dry area at room temperature.

**Storage:** Upon receipt, store the kit at 2 – 8°C. DO NOT FREEZE. Under these conditions the reagents will retain their activity until the date shown on the bottle labels. After use, return the kit to the refrigerator, storing bottles in an upright position.

## **SPECIMEN COLLECTION AND PREPARATION**

Staphylococcal isolates grown on nonselective media, such as blood agar, should be used.<sup>3</sup> Selective media, such as Baird-Parker or Mannitol Salt Agar, may also give satisfactory results, but interpretation is more difficult due to weak positive reactions and the tendency to produce stringy reactions. Subculture of isolates onto nonselective media is recommended. Test only fresh cultures or subcultures (18 – 36 h incubation).

## **PROCEDURES**

Review “Precautions” and “Specimen Collection and Preparation” prior to performing procedures. The testing area, reagents, test specimens and test components should be at room temperature when used.

**Materials Provided:** All materials listed under “Reagents”, and test disposables and accessories.

**Materials Required But Not Provided:** Timer, microbiological loop or wooden mixing sticks, Bunsen burner, disinfectant, Gram stain reagents and quality control organisms.

Also required are the necessary equipment and labware used for preparation, storage and handling of specimens.

**Test Procedure:** If fewer than six tests are performed, the test card may be cut with scissors and the unused portions saved for later use.

1. Mix the latex reagent by shaking; expel any latex from the dropper for complete mixing.
2. Dispense 1 drop of Test Latex onto one of the circles on the reaction card and 1 drop of Control Latex onto another circle.
3. Using a microbiological loop pick up and smear 5 suspect colonies onto the Test Latex-containing circle and mix this into the Test Latex reagent. Spread to cover the circle.
4. Repeat step 3 for the Control Latex.
5. Pick up and hand rock the card for up to 20 sec and observe for agglutination under normal lighting conditions. Read macroscopically; do not use a magnifying glass.
6. Dispose of the reaction card in an appropriate biohazard container.
7. Re-cap the bottles and return to the refrigerator.

**User Quality Control:** Perform Daily Quality Control each day the kit is used. Follow the methods given in “Test Procedure”.

1. *Positive Control:* Use a known *S. aureus* strain such as ATCC™ 25923. Ensure that agglutination occurs within 20 sec.

2. **Negative Control:** Use a known *S. epidermidis* strain such as ATCC 12228. Ensure that the reagent remains smooth and unagglutinated for the entire 20 sec of the test.

**NOTE:** Do not use the reagents if reactions with the control organisms are incorrect.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

## INTERPRETATION OF TEST RESULTS

**Positive Results:** A positive result is obtained if agglutination of the blue latex particles occurs within 20 sec in the test circle, with no agglutination in the control circle. The result is positive when there is noticeable clearing of the blue background in the test latex. This indicates the presence of *S. aureus*.

**Negative Results:** A negative result is obtained if no agglutination occurs and a smooth suspension remains at 20 sec in the test circle. The result is negative when there is no noticeable clearing of the blue background in the test latex.

Reactions occurring after 20 sec should be ignored.

Occasional granular or stringy reactions may be seen due to the particular nature of the test. If granular or stringy reactions occur, they should be interpreted using the following criteria:

The result is **POSITIVE** when there is a noticeable clearing of the blue background.

The result is **NEGATIVE** when there is *no* noticeable clearing of the blue background.

**Uninterpretable Results:** The test is uninterpretable if the control reagent shows agglutination or autoagglutination.

## LIMITATIONS OF THE PROCEDURE

Some staphylococci other than *S. aureus* may give positive coagulase results. These strains include *S. hyicus*, *S. intermedius*, *S. lugdunensis*, and *S. schleiferi*, which may also react in rapid latex tests. If necessary, these species would require identification using biochemical test procedures. Both *S. intermedius* and *S. hyicus* are rarely isolated from human specimens.<sup>1</sup>

The amount of protein A expressed by *S. aureus* is dependent on the medium, and other growth conditions. Rare strains of *S. aureus* that do not produce coagulase but still express some protein A on their cell surface may yield a weak positive reaction with the **BBL Staphyloslide** Latex Test. If such a strain is suspected, alternative identification methods should be used.

Studies have shown that a limited number of MRSA strains may produce weak levels of clumping factor and protein A. These strains may fail to react in the latex test.<sup>4</sup>

Staphylococci isolated from urine specimens which give a weak positive<sup>5</sup> or stringy result with the **BBL Staphyloslide** Latex Test may be *Staphylococcus*

*saprophyticus*. Further identification of such isolates may be conducted using biochemical tests and novobiocin sensitivity (*S. saprophyticus* is resistant to novobiocin).

Autoagglutination, where there is a partial clearing of the blue background, may occur if rough strains of staphylococci are tested or if a culture is incubated beyond 36 h. If the test colonies appear rough or autoagglutination is suspected, the organism should be subcultured and retested at 18 – 24 h. If the organism continues to produce autoagglutination, the test result is uninterpretable.

Some streptococci and possibly other organisms possessing immunoglobulin or plasma binding factors may react in the latex test, and some species such as *Escherichia coli* are able to nonspecifically agglutinate latex particles.<sup>6,7</sup> To overcome these potential nonspecific results, a Gram stain should be performed so only typical staphylococci are tested.

The use of high salt media may result in a weaker reaction. Colonies taken from these media may be more difficult to emulsify and may give rise to slightly stringy reactions.

### PERFORMANCE CHARACTERISTICS

In three clinical comparisons of the performance of the BBL Staphyloslide Latex Test, the sensitivity of the Latex Test was 100% (530/530) and the specificity was 99% (379/383).

	BBL Staphyloslide Latex		Commercial Latex Test		Tube Coagulase		No. of Strains Tested
	+	-	+	-	+	-	
<i>S. aureus</i> <sup>(a)</sup>	530	0	530	0	526	4 <sup>(c)</sup>	530
Non- <i>S. aureus</i>	4 <sup>(b)</sup>	379	4 <sup>(b)</sup>	379	0	383	383

- To be classified as *S. aureus*, a positive result must be obtained in two or more established methods; e.g., slide coagulase, tube coagulase, DNase or biochemical tests.
- All strains were identified biochemically as *S. saprophyticus*.
- All discrepant strains were identified as *S. aureus* as determined by two or more of the following tests: slide coagulase, red cell agglutination, latex agglutination or biochemical methods.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and Staphyloslide are trademarks of Becton, Dickinson and Company © 2010 BD.
Benex Limited