## POSTUPCI KONTROLE KVALITETE

## I UVOD

**BD BBL *Clostridium difficile* Selective Agar (CDSA)** (Selektivni agar za *Clostridium difficile* (CDSA)) preporučuje se kao selektivni i diferencijalni medij za primarnu izolaciju bakterije *Clostridium difficile* iz fekalnih uzoraka.

## II POSTUPAK ISPITIVANJA UČINKOVITOSTI

1. Stavite pločice na sobnu temperaturu na 6–24 h prije upotrebe u anaerobnom sustavu **BD GasPak EZ**.
2. Inokulirajte reprezentativne uzorke s dolje navedenim kulturama.
  - a. Za obligatne anaerobe razmažite inokulat 1 µL (0,001 mL) iz 48–72 h kulture nasjeckanog mesnog bujona (Chopped Meat Broth) koji je bio inkubiran 2 dana na temperaturi 35–37 °C.
  - b. Za *E. coli* i *P. mirabilis* razmažite inokulat 1 µL (0,001 mL) iz 4–5 h kulture **BD BBL Trypticase** sojinog bujona razrijeđenog da daje 10<sup>6</sup>–10<sup>7</sup> CFU/mL.
  - c. Inkubirajte pločice pri temperaturi od 36 °C ± 1 °C u anaerobnoj atmosferi.
  - d. Uključite pločice prethodno testirane serije **BD Trypticase** sojinog agara s 5% ovčje krvi kao kontrole za inhibirane sojeve.
3. Nakon 48–72 h na pločicama pregledajte rast, boju kolonije, fluorescenciju i selektivnost.
4. Očekivani rezultati

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222228 **BD BBL *Clostridium difficile* Selective Agar**

## XIV REFERENCE

1. Bartlett, J.G., T.W. Chang, M. Gurwith, S.L. Gorbach, and A.B. Onderdonk. 1978. Antibiotic-associated pseudomembranous colitis due to toxin-producing clostridia. *N. Engl. J. Med.* 298:531-534.
2. Koransky, J.R., S.D. Allen, and V.R. Dowell, Jr. 1978. Use of ethanol for isolation of spore forming microorganisms. *Appl. Environ. Microbiol.* 35:762-765.
3. Kafiz, S., and C.L. Oakley. 1976. *Clostridium difficile*: isolation and characteristics. *J. Med. Microbiol.* 9:129-137.
4. Willey, S.H., and J.G. Bartlett. 1979. Cultures for *Clostridium difficile* in stools containing a cytotoxin neutralized by *Clostridium sordellii* antitoxin. *J. Clin. Microbiol.* 10:880-884.
5. George, W.L., V.L. Sutter, D. Citron, and S.M. Finegold. 1979. Selective and differential medium for isolation of *Clostridium difficile*. *J. Clin. Microbiol.* 9:214-219.
6. Wust, J., N.M. Sullivan, U. Hardegger, and T.D. Wilkins. 1982. Investigation of an outbreak of antibiotic-associated colitis by various typing methods. *J. Clin. Microbiol.* 16:1096-1101.
7. Dowell, V.R., Jr. 1981. Media for the selective isolation of *Clostridium difficile*. *Les Anaerobes Microbiologie-Pathologie. Symposium International.* New York.
8. Isenberg, H.D., F.D. Schoenkecht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
9. Miller, J.M., and H.T. Holmes. 1995. Specimen collection, transport, and storage, p. 19-32. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
10. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
11. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
12. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
13. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
14. Allen, S.D., J.A. Siders, and L.M. Marler. 1985. Isolation and examination of anaerobic bacteria, p. 413-433. In E.H. Lennette, A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.), *Manual of clinical microbiology*, 4th ed. American Society for Microbiology, Washington, D.C.
15. Martin, W.J. 1971. Practical method for isolation of anaerobic bacteria in the clinical laboratory. *Appl. Microbiol.* 22:1168-1171.
16. George, W.L., R.D. Rolfe, and S.M. Finegold. 1982. *Clostridium difficile* and its cytotoxin in feces of patients with antimicrobial-associated diarrhea and miscellaneous conditions. *J. Clin. Microbiol.* 15:1049-1053.
17. Baron, E.J., Assessment of currently available laboratory tests for *Clostridium difficile*-associated diarrhea. *Clin. Microbiol. Newsl.* 11:118-120.
18. Dowell, V.R., Jr., and T.M. Hawkins. 1987. Laboratory methods in anaerobic bacteriology. CDC laboratory manual. HHS Publication No. (CDC) 87-8272. Centers for Disease Control, Atlanta.
19. Holdeman, L.V., E.P. Cato, and W.E.C. Moore (ed.). 1977. Anaerobe laboratory manual, 4th ed. Virginia Polytechnic Institute and State University, Blacksburg.
20. Rodloff, A.C., P.C. Applebaum, and R.J. Zabransky. 1991. Cumitech 5A, Practical anaerobic bacteriology. Coordinating ed., A.C. Rodloff, American Society for Microbiology, Washington, D.C.
21. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.
22. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. 1994. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
23. Allen, S.D., C.L. Emery, and J.A. Siders. 1999. *Clostridium*, p. 654-671. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), *Manual of clinical microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.

Tehnički servis i podrška: obratite se lokalnom predstavniku tvrtke BD ili posjetite [www.bd.com](http://www.bd.com).

Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.