POSTUPCI KONTROLE KVALITETE**I UVOD**

Kligler Iron Agar (Agar željeza Kligler) pomaže u diferencijaciji gram negativnih enterobacila na temelju njihove sposobnosti da fermentiraju dekstrozu i laktozu i stvaraju sulfide.

II POSTUPAK ISPITIVANJA UČINKOVITOSTI

1. Inokulirajte reprezentativne uzorke s dolje navedenim kulturama.
 - a. Pomoću kultura kosog sojinog agara **Trypticase** starih 18 – 24 h inokulirajte epruvete inokulacijskom iglom tako da ubodete kraj epruvete i razmažete cijelu površinu kosog agara.
 - b. Inkubirajte epruvete s malo odvjanim čepovima na 35 ± 2 °C u aerobnoj atmosferi.
2. Pregledajte epruvete nakon 18 – 24 h i potražite rast i reakcije.
3. Očekivani rezultati

Organizmi	ATCC	Kosi agar	Kraj	Plin	H ₂ S
* <i>Escherichia coli</i>	25922	Kiseo	Kiseo	+	–
* <i>Salmonella enterica</i> podvrst. <i>enterica</i> serotip Typhimurium	14028	Lužnat	Kiseo	+/-	+
* <i>Shigella flexneri</i>	12022	Lužnat	Kiseo	–	–
<i>Pseudomonas aeruginosa</i>	27853	Lužnat	Lužnat	–	–

*Preporučeni soj organizama za korisničku kontrolu kvalitete.

III DODATNA KONTROLE KVALITETE

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XIV REFERENCE

- Russell, F.F. 1911. The isolation of typhoid bacilli from urine and feces with the description of a new double sugar tube medium. *J. Med. Res.* 25:217-229.
- Kligler, I.J. 1917. A simple medium for the differentiation of members of the typhoid-paratyphoid group. *Am. J. Public Health.* 7:1042-1044.
- Kligler, I.J. 1918. Modifications of culture media used in the isolation and differentiation of typhoid, dysentery and allied bacilli. *J. Exp. Med.* 28:319-322.
- Bailey, S.F., and N.I. Lacy. 1927. A modification of the Kligler lead acetate medium. *J. Bacteriol.* 13:183-189.
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.). 2003. *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. *Bailey & Scott's diagnostic microbiology*, 11th ed. Mosby, Inc., St. Louis.
- Ewing, W.H. 1986. *Edwards and Ewing's identification of the Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc. New York.
- MacFaddin, J.F. 1985. *Media for isolation-cultivation-identification-maintenance of medical bacteria*, vol. I. Williams & Wilkins, Baltimore.
- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Year Book, Inc., St. Louis.
- Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. *Bergey's Manual™ of determinative bacteriology*, 9th ed. Williams & Wilkins, Baltimore.

Tehnički servis i podrška na BD Diagnostics: obratite se lokalnom predstavniku tvrtke BD ili posjetite www.bd.com/ds.

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

 Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL, GasPak and Trypticase are trademarks of Becton, Dickinson and Company. ©2014 BD.