BD MacConkey II Agar

NAMJENA

BD MacConkey II Agar (Agar MacConkey II) selektivna je i diferencijalna podloga za izolaciju i diferencijaciju *Enterobacteriaceae* i mnoštvo drugih rodova gram negativnih bakterija iz kliničkih uzoraka.

NAČELA I OBJAŠNENJE POSTUPKA

Mikrobiološka metoda.

U današnje vrijeme mnogo je hranjivih podloga dostupno za izolaciju, uzgoj i identifikaciju *Enterobacteriaceae* i određenih nefermentirajućih bakterija. Jednu od najranijih takvih podloga razvio je MacConkey i objavio 1900 i 1905. godine.^{1,2} Ta formula temelji se na poznavanju činjenice da se žučne soli talože uslijed djelovanja kiselina, a da određeni enterični mikroorganizmi fermentiraju laktozu, dok drugi nemaju tu sposobnost. Kasnije je podloga mijenjana nekoliko puta.^{3,4}

Agar MacConkey samo je blago selektivan budući da je koncentracija žučnih soli koje inhibiraju gram pozitivne mikroorganizme niska u usporedbi s drugim enteričnim pločastim podlogama. Ova podloga preporučuje se za upotrebu s kliničkim uzorcima za koje se pretpostavlja da sadrže miješanu mikrobnu floru, kao što su urin, respiratorni uzorci, rane i drugo, zbog toga što omogućuje preliminarno grupiranje enteričnih i drugih gram negativnih bakterija u one koje fermentiraju laktozu i u one koje ju ne fermentiraju.^{5,6} Agar MacConkey također se koristi i za mikrobiološko ispitivanje hrane.⁷

Formula agara MacConkey II osmišljena je da poboljša inhibiciju množenja vrsta *Proteus* radi dobivanja konačne diferencijacije na one bakterije koje fermentiraju laktozu i one koje je ne fermentiraju i za izrazit rast enterobakterija.

Peptoni snabdjevaju **BD MacConkey II Agar** hranjivim tvarima. Kristalno ljubičasta boja inhibira gram pozitivne bakterije, osobito enterokoke i stafilokoke. Diferencijacija enteričnih mikroorganizama postiže se kombinacijom laktoze i neutralnog crvenog pH indikatora. Stvaraju se bezbojne ili ružičaste do crvene kolonije ovisno o sposobnosti izolata da fermentira ugljikohidrate.

REAGENSI

BD MacConkey II Agar

Formula* po litri pročišćene vode

Pankreatična digestija želatine	17,0 g
Pankreatična digestija kazeina	1,5
Peptična digestija životinjskog tkiva	1,5
Laktoza	10,0
Žučne soli	1,5
Natrijev klorid	5,0
Neutralna crvena boja	0,03
Kristalno ljubičasta boja	0,001
Agar	13,5

pH 7,1 ± 0,2

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REFERENCE

1. MacConkey, A.T. 1900. Note on a new medium for the growth and differentiation of the *Bacillus coli communis* and the *Bacillus typhi abdominalis*. The Lancet, Part II:20.
2. MacConkey, A. 1905. Lactose-fermenting bacteria in faeces. J. Hyg. 5:333-379.
3. Levine, M., and H.W. Schoenlein. 1930. A compilation of culture media for the cultivation of microorganisms. The Williams & Wilkins Company, Baltimore.
4. MacFaddin, J.F. 1985. Media for isolation-cultivation- identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
5. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
6. Farmer III, J.J. 2003. *Enterobacteriaceae*: introduction and identification. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
7. Downes, F.P., and K. Ito. 2001. Compendium of methods for the microbiological examination of foods. 4th edition. American Public Health Association (APHA). Washington, D.C. USA.
8. Thomson, R.B., and J.M. Miller. 2003. Specimen collection, transport, and processing: bacteriology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
9. Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
10. Mazura-Reetz, G., T.R. Neblett, and J.M. Galperin. 1979. MacConkey agar: CO₂ vs. ambient incubation, abstr. C 179, p. 339. Abstr. 79th Annu. Meet. Am. Soc. Microbiol. 1979.

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.

© 2014 BD

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Becton Dickinson GmbH