**BD Sabouraud Glucose Agar •**  
**BD Sabouraud Agar with Chloramphenicol •**  
**BD Sabouraud Agar with Gentamicin and Chloramphenicol •**  
**BD Sabouraud Agar with Penicillin and Streptomycin •**

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Priloženi materijal

**BD Sabouraud Glucose Agar,**

**BD Sabouraud Agar with Chloramphenicol,**

**BD Sabouraud Agar with Gentamicin and Chloramphenicol ili**

**BD Sabouraud Agar with Penicillin and Streptomycin,** svi isporučeni u pločicama **Stacker** od 90 mm. Mikrobiološki kontrolirano.

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## REFERENCE

1. Sabouraud, R. 1892. Contribution a l'etude de la trichophytie humaine. Etude clinique, microscopique et bacteriologique sur la pluralité des trichophyton de l'homme. Ann. Dermatol. Syphil. 3: 1061-1087.
2. Haley, L.D., J. Trandel, and M.B. Coyle. 1980. Cumitech 11, Practical methods for culture and identification of fungi in the clinical microbiology laboratory. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.
3. Ajello, L., L.K. Georg, W. Kaplan, and L. Kaufman. 1963. CDC laboratory manual for medical mycology. PHS Publication No. 994, U.S. Government Printing Office, Washington, D.C.
4. MacFaddin, J.F. 1985. Media for isolation-cultivation- identification-maintenance of medical bacteria. vol. I. Williams & Wilkins, Baltimore.
5. Sutton, D.A. 2003. Specimen Collection, transport, and processing: mycology. *In*: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
6. U.S. Pharmacopeial Convention, Inc. 2009. The U.S. Pharmacopeia 32/The national formulary 27--2009. U.S. Pharmacopeial Convention, Inc., Rockville, Md. USA
7. Larone, D.H. 2002. Medically important fungi: a guide to identification. 4<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
8. Summerbell, R.C. 2003. *Trichophyton, Microsporum, Epidermophyton*, and agents of superficial mycoses. *In*: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
9. Kwon-Chung, K.J., and J.E. Bennett. 1992. Medical mycology. Lea & Febiger, Philadelphia.
10. Pfaller, M.A., and R.A. Tenover (section ed.). 2003. Mycology. *In*: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
11. Brun, S., et al. 2001. Evaluation of five commercial Sabouraud gentamicin-chloramphenicol agar media. Eur. J. Clin. Microbiol. Infect. Dis. 20: 718-723.
12. Council of Europe, 2008. European Pharmacopoeia, 6.1. European Pharmacopoeia Secretariat. Strasbourg/France.

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### BD Sabouraud Agar with Chloramphenicol

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### BD Sabouraud Agar with Gentamicin and Chloramphenicol

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### BD Sabouraud Agar with Penicillin and Streptomycin

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection  
© 2019 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

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### Becton Dickinson GmbH