BD Haemophilus Test Medium Agar (HTM)

NAMJENA

BD Haemophilus Test Medium Agar (Agar s testnim medijem Haemophilus) (HTM) koristi se u postupku protumikrobnog disk difuzijskog ispitivanja osjetljivosti za kliničke izolate vrste *Haemophilus* kao što je opisano u odobrenom standardu M2 koji je objavio Institut za kliničke i laboratorijske standarde (CLSI).¹

NAČELA I OBJAŠNJENJE POSTUPKA

Mikrobiološka metoda.

Agar Mueller Hinton s dodatkom 1% hemoglobina i obogaćenja 1% IsoVitaleX (čokoladni agar Mueller Hinton) podloga je koja se ranije preporučivala za *Haemophilus influenzae*.² U opsežnoj studiji koju su proveli Jorgensen i kolege razvijen je agar s testnim medijem Haemophilus (Haemophilus Test Medium HTM).^{3,4} Ta podloga je agar Mueller Hinton s faktorom (hemin ili hematin), faktorom V (nikotinamid adenin dinukleotid, NAD) i ekstraktom kvasca.

Najveća prednost agara HTM u usporedbi s čokoladnim agarom Mueller Hinton optička je čistoća koja omogućuje mjerenja promjera zona s dna posudice kakav i je standardni postupak ispitivanja za nezahtjevne organizme na agaru Mueller Hinton. Osim toga, agar HTM sadrži nisku razinu timidina i zbog toga je prikladan za ispitivanje trimetoprima/sulfametoksazola.

Kriteriji za interpretaciju protumikrobnog ispitivanja osjetljivosti navedeni su u CLSI-jevom dokumentu M100 (M2) koji je uključen u CLSI-jev dokument M2.¹ Taj dokument treba konzultirati za dodatne pojedinosti.

BD Haemophilus Test Medium Agar (Haemophilus Test Medium) sastoji se od agara Mueller Hinton u koji je dodan faktor X (hemin ili hematin), faktor V (nikotinamid adenin dinukleotid, NAD) i ekstrakt kvasca.

REAGENSI

BD Haemophilus Test Medium Agar (HTM)

Formula* po litri pročišćene vode

Goveđi ekstrakt	2,0 g
Kiseli hidrolizat kazeina	17,5
Škrob	1,5
Agar	17,0
Ekstrakt kvasca	5,0
Hematin	0,015
Nikotinamid-adenin-dinukleotid (NAD)	0,015

pH 7,3 ± 0,1

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REFERENCE

1. Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS). Approved standard: M2. Performance standards for antimicrobial disk susceptibility tests. CLSI, Wayne, PA, USA. Search for latest version at www.clsi.org
2. Thornsberry, C., T.L. Gavan, E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C
3. Jorgensen, J.H., J.S. Redding, L.A. Maher, and A.W. Howell. 1987. Improved medium for antimicrobial susceptibility testing of *Haemophilus influenzae*. J. Clin. Microbiol. 25:2105-2113.
4. Jorgensen, J.H., A.W. Howell, and L.A. Maher. 1990. Antimicrobial susceptibility testing of less commonly isolated *Haemophilus* species using *Haemophilus* test medium. J. Clin. Microbiol. 28:985-988.
5. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. J. Clin. Microbiol. 17:450-457.
6. Marsik, F., G. Evans, J. Fowler, and L. Thompson. 1989. Comparison of the BBL Prompt system, Abbott A-Just, and visual method for the preparation of *Haemophilus influenzae* inoculum for the Bauer-Kirby procedure, abstr. C-67, p. 404. Abstr. 89th Annu. Meet. Am. Soc. Microbiol. 1989.
7. Doern, G.V., J.H. Jorgensen, C. Thornsberry, D.A. Preston, T.A. Tubert, J.S. Redding, and L.A. Maher. 1988. A national collaborative study of the prevalence of antimicrobial resistance among clinical isolates of *Haemophilus influenzae*. Antimicrob. Agents Chemother. 32:180-185.
8. Službeni podaci. BD Diagnostic Systems, Sparks, MD, SAD.
9. National Committee for Clinical Laboratory Standards. 1990. Approved standard: M2-A4. Performance standards for antimicrobial disk susceptibility tests, 4th ed. National Committee for Clinical Laboratory Standards, Villanova, Pa.
10. Azemun, P., T. Stull, M. Roberts, and A.L. Smith. 1981. Rapid detection of chloramphenicol resistance in *Haemophilus influenzae*. Antimicrob. Agents Chemother. 20:168-170.
11. Doern, G.V., G.S. Daum, and T.A. Tubert. 1987. In vitro chloramphenicol susceptibility testing of *Haemophilus influenzae*: disk diffusion procedures and assays for chloramphenicol acetyltransferase. J. Clin. Microbiol. 25:1453-1455.
12. National Committee for Clinical Laboratory Standards. 1984. Approved standard: M2-A3. Performance standards for antimicrobial disk susceptibility tests, 3rd ed. National Committee for Clinical Laboratory Standards, Villanova, Pa.
13. National Committee for Clinical Laboratory Standards. 1987. Second Informational Supplement: M100-S2. Performance standards for antimicrobial susceptibility testing. National Committee for Clinical Laboratory Standards, Villanova, Pa.
14. Doern, G.V. and R.N. Jones. 1988. Antimicrobial susceptibility testing of *Haemophilus influenzae*, *Branhamella catarrhalis*, and *Neisseria gonorrhoeae*. Antimicrob. Agents Chemother. 32:1747-1753.
15. Woods, G.L., and J.A. Washington. 1995. Antibacterial susceptibility tests: dilution and disk diffusion methods, p. 1327-1341. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C. USA.

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<http://www.bd.com/europe/regulatory/>

Prompt is a trademark of 3M

ATCC is a trademark of the American Type Culture Collection

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2013 BD

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Becton Dickinson GmbH