## **BD Mueller Hinton II Agar** **BD Mueller Hinton II Agar 150 mm** **BD Mueller Hinton II Agar, Square**

### **NAMJENA**

**BD Mueller Hinton II Agar** (Agar BD Mueller Hinton II), dostupan u nekoliko formata pločica, koristi se u standardiziranom disk difuzijskom postupku za utvrđivanje osjetljivosti brzorastućih aerobnih organizama na protumikrobnim agensima prema standardima Instituta za kliničke i laboratorijske studije (CLSI) i Europskog povjerenstva za ispitivanje antimikrobne osjetljivosti European Committee on Antimicrobial Susceptibility Testing (EUCAST).<sup>1, 2</sup>

### **NAČELA I OBJAŠNJENJE POSTUPKA**

Mikrobiološka metoda.

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## REFERENCE

1. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
2. Matuschek, E., D.F. Brown, and G. Kahlmeter. Development of the EUCAST disk diffusion antimicrobial susceptibility testing method and its implementation in routine microbiology laboratories. *Clin. Microbiol. Infect.* 2014; 20 (4): 255-66.
3. Ryan, K.J., F.D. Schoenknecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice* 5:91-100.
4. Barry, A.L., F. Garcia, and L.D. Thrupp. 1970. An improved single-disk method for testing the antibiotic susceptibility of rapidly-growing pathogens. *Am. J. Clin. Pathol.* 53:149-158.
5. Mueller, J.H., and J. Hinton. 1941. A protein-free medium for primary isolation of the gonococcus and meningococcus. *Proc. Soc. Exp. Biol. Med.* 48:330-333.
6. Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta Pathol. Microbiol. Scand. Sec. B, Suppl.* 217.
7. CLSI. Approved standard: M02-A12 - Performance standards for antimicrobial disk susceptibility tests. CLSI, Wayne, PA, USA. *Search for latest version at [www.clsi.org](http://www.clsi.org).*
8. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing. *Search for latest version at <http://www.eucast.org>.*
9. Woods, G.L., and J.A. Washington. 1995. Antibacterial susceptibility tests: dilution and disk diffusion methods, p. 1327-1341. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.C. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, DC.
10. CLSI - Performance Standards for Antimicrobial Susceptibility. CLSI supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute. *Search for latest version at [www.clsi.org](http://www.clsi.org).*
11. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. *Search for latest version at <http://www.eucast.org>.*
12. Thornsberry, C., T.L. Gavan, and E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society of Microbiology, Washington, DC.

13. Koch, A.E., and J.J. Burchall. 1971. Reversal of the antimicrobial activity of trimethoprim by thymidine in commercially prepared media. *Appl. Microbiol.* 22:812-817.
14. Ferone, R., S.R.M. Bushby, J.J. Burchall, W.D. Moore, and D. Smith. 1975. Identification of Harper-Cawston factor as thymidine phosphorylase and removal from media of substances interfering with susceptibility testing to sulfonamides and diaminopyrimidines. *Antimicrob. Agents Chemother.* 7:91-98.
15. Reller, L.G., F.D. Schoenknecht, M.A. Kenny, and J.C. Sherris. 1974. Antibiotic susceptibility testing of *Pseudomonas aeruginosa*: selection of a control strain and criteria for magnesium and calcium content in media. *J. Infect. Dis.* 130:454-463.
16. Pollock, H.M., B.H. Minshew, M.A. Kenny, and F.D. Schoenknecht. 1978. Effect of different lots of Mueller-Hinton Agar on the interpretation of the gentamicin susceptibility of *Pseudomonas aeruginosa*. *Antimicrob. Agents Chemother.* 14:360-367.
17. D'Amato, R.F., and C. Thornsberry. 1979. Calcium and magnesium in Mueller-Hinton agar and their influence on disk diffusion susceptibility results. *Current Microbiol.* 2:135-138.
18. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. *Search for latest version at <http://www.eucast.org>.*
19. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 6.0, 2016. *<http://www.eucast.org>.*
20. CLSI. Performance Standards for Antimicrobial Susceptibility Testing. 26<sup>th</sup> ed. CLSI supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2016.
21. Wegner, D.L., C.R. Mathis, and T.R. Neblett. 1976. Direct method to determine the antibiotic susceptibility of rapidly growing blood pathogens. *Antimicrob. Agents Chemother.* 9:861-862.
22. Johnson, J.E., and J.A. Washington II. 1976. Comparison of direct and standardized antimicrobial susceptibility testing of positive blood cultures. *Antimicrob. Agents Chemother.* 10:211-214.
23. Waterworth, P.M., and M. Del Piano. 1976. Dependability of sensitivity tests in primary culture. *J. Clin. Pathol.* 29:179-184.
24. Hollick, G.E., and J.A. Washington II. 1976. Comparison of direct and standardized disk diffusion susceptibility testing of urine cultures. *Antimicrob. Agents Chemother.* 9:804-809.
25. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
26. Hindler, J.A., and C.B. Anderbied. 1985. Effect of the source of Mueller-Hinton agar and resistance frequency on the detection of methicillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.* 21:205-210.
27. Maskell, R., O.A. Okubadejo, R.H. Payne, and L. Pead. 1977. Human infections with thymine-requiring bacteria. *J. Med. Microbiol.* 11:33-45.
28. Haltiner, R.C., P.C. Migneault, and R.G. Robertson. 1980. Incidence of thymidine-dependent enterococci detected on Mueller-Hinton agar with low thymidine content. *Antimicrob. Agents Chemother.* 18:365-368.
29. Murray, B.E. 1990. The life and times of the *Enterococcus*. *Clin. Microbiol. Rev.* 3:46-65.
30. CLSI. Approved standard: M7. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically. CLSI, Wayne, PA, USA. *Search for latest version at [www.clsi.org](http://www.clsi.org).*
31. Jorgensen, J.H., and J.D. Turnidge. 2003. Susceptibility test methods: dilution and disk diffusion methods. *In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Tenover, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.*
32. Washington, J.A., and G.L. Woods. 1995. Antimicrobial susceptibility tests: dilution and disk diffusion methods. P. 1327-1341. *In Muarry, P.R., Baron, E.J., Tenover, M.A., Tenover, F.C., and Tenover, R.H. (ed.), Manual of clinical microbiology, 6<sup>th</sup> ed. American Society for Microbiology, Washington, D.C. 1995.*

33. Neumann, M.A., D.F. Sahn, C. Thornsberry, J.E. McGowan, Jr. Cumitech 6A, New developments in antimicrobial agent susceptibility testing: a practical guide. Coordinating ed., J.E. McGowan, Jr. American Society of Microbiology, Washington, D.C.1991.
34. Data on file. Becton Dickinson.
35. Skov, R., A.R. Larsen, A. Kearns, M. Holmes, C. Teale, G. Edwards, and R. Hill. Phenotypic detection of *mecC*-MRSA: cefoxitin is more reliable than oxacillin. 2014. J. Antimicrob. Chemother. 69:133-135.

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© 2017. BD. BD, logotip BD i svi drugi zaštitni znakovi u vlasništvu su tvrtke Becton, Dickinson and Company.
### **Becton Dickinson GmbH**