## **BD MacConkey II Agar / Columbia CNA Agar Improved II with 5% Sheep Blood (Biplate)**

### **NAMJENA**

**BD MacConkey II Agar / Columbia CNA Agar Improved II with 5% Sheep Blood (Biplate)** (agar BD MacConkey II / poboljšani agar Columbia CNA II s 5% ovčje krvi (dvostruka ploča)) poboljšana je podloga koja se koristi za selektivnu izolaciju iz gram-pozitivnih i gram-negativnih bakterija iz kliničkih uzoraka.

### **NAČELA I OBJAŠNJENJE POSTUPKA**

Mikrobiološka metoda.

MacConkey Agar jedna je od najranijih formula podloge (objavio ju je 1900. MacConkey) za izolaciju, uzgoj i identifikaciju *Enterobacteriaceae* i određenih nefermentirajućih bakterija. Ta je podloga naknadno izmijenjena nekoliko puta.<sup>1,2</sup>

Formula za MacConkey II Agar osmišljena je 1987. da poboljša inhibiciju množenja vrsta *Proteus* radi dobivanja konačne diferencijacije na one bakterije koje fermentiraju laktozu i one koje je ne fermentiraju i za izrazit rast enterobakterija. Peptoni snabdijevaju agar MacConkey II Agar hranjivim tvarima. Kristalno ljubičasta boja služi inhibiciji gram-pozitivnih bakterija, osobito enterokoka i stafilokoka. Diferencijacija enteričnih mikroorganizama postiže se kombinacijom laktoze i neutralnog crvenog pH indikatora. Stvaraju se bezbojne ili ružičaste do crvene kolonije, ovisno o sposobnosti izolata da fermentira ugljikohidrate.<sup>3-5</sup>

Godine 1966. Ellner i dr. prijavili su razvoj nove formule krvnog agara koji je nazvan Columbia Agar<sup>6</sup> Ova podloga na kojoj se dobivaju veće kolonije i bujniji rast nego na usporednim bazama krvnih agara, koristi se za podloge koje sadrže krv i za selektivne formule. Ellner i dr. otkrili su da podloga koja sadrži 10 mg kolistina i 15 mg nalidiksične kiseline po litri u agarnoj bazi Columbia, obogaćena s 5% ovčje krvi, potiče rast stafilokoka, hemolitičkih streptokoka i enterokoka, a inhibira rast vrsta *Proteus*, *Klebsiella* i *Pseudomonas*.<sup>6</sup> Tijekom godina, povećala se otpornost bakterija na protumikrobne agense. To je posebno svojstveno za gram-negativne štapiće koji se trebaju inhibirati, no često omogućuju rast na agaru Columbia CNA Agar with 5% Sheep Blood. U agaru Columbia CNA Agar Improved II with 5% Sheep Blood, mala količina aztreonama održava dobru selektivnost podloge, a koncentracija nalidiksične kiseline smanjena je na 5,5 mg/L kako bi se povećala izolacija gram-pozitivnih koka, posebno stafilokoka. Koncentracija kolistina ostala je ista. Aztreonam je monobaktam koji djeluje samo protiv većine gram-negativnih bakterija, dok na gram-pozitivne organizme ne djeluje.<sup>7-9</sup> Ovčja krv omogućuje otkrivanje hemolitičkih reakcija koje su osobito važne u vjerojatnoj dijagnozi streptokoka.<sup>10</sup>

Glavna prednost agara Columbia CNA Agar Improved II with 5% Sheep Blood u odnosu na agar Columbia CNA Agar with 5% Sheep Blood poboljšani je rast stafilokoka koji se češće otkrivaju nakon 18 do 24 sata inkubacije te bolja inhibicija gram-negativnih bakterija, posebno *Proteus* spp.

Kombinacija ove dvije podloge na dvostrukoj ploči (**BD MacConkey II Agar / agar Columbia CNA Agar Improved II with 5% Sheep Blood**) koristi se za selektivnu izolaciju iz gram-pozitivnih i gram-negativnih bakterija iz kliničkih uzoraka.

## REAGENSI

### BD MacConkey II Agar / Columbia CNA Agar Improved II with 5% Sheep Blood (Biplate)

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## REFERENCE

1. Levine, M., and H.W. Schoenlein. 1930. A compilation of culture media for the cultivation of microorganisms. The Williams & Wilkins Company, Baltimore.
2. MacFaddin, J.F. 1985. Media for isolation-cultivation- identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
3. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
4. Farmer, J.J., III., K.D. Boatwright, and J.M. Janda 2007. *Enterobacteriaceae*: introduction and identification. In: Murray, P. R., E. J. Baron, M.L. Landry, J.H. Jorgensen, and M. A. Pfaller (ed.). Manual of clinical microbiology, 9<sup>th</sup>ed. American Society for Microbiology, Washington, D.C.
5. Downes, F.P., and K. Ito. 2001. Compendium of methods for the microbiological examination of foods. 4<sup>th</sup> edition. American Public Health Association (APHA). Washington, D.C. USA.
6. Ellner, P.D., C.J. Stoessel, E. Drakeford, and F. Vasi. 1966. A new culture medium for medical bacteriology. Am. J. Clin. Pathol. 45: 502-504.
7. Wood, W., G. Harvey, E.S. Olson, and T.M. Reid. 1993. Aztreonam selective agar for Gram positive bacteria. J. Clin. Pathol. 46: 769-771.
8. Wiedemann, B., and B. A. Atkinson. 1986. Susceptibility to antibiotics: species incidence and trends. In: Lorian, V. (ed.), Antibiotics in Laboratory medicine, p. 962-1208. Williams and Wilkins, Baltimore, USA.
9. von Graevenitz, A. 1986. Use of antimicrobial agents as tools in epidemiology, identification, and selection of microorganisms. In: Lorian, V. (ed.), Antibiotics in Laboratory medicine, p. 723-738. Williams and Wilkins, Baltimore, USA.
10. Spellerberg, B., Brandt, C. 2007. *Streptococcus*. In: Murray, P. R., E. J. Baron, M.L. Landry, J.H. Jorgensen, and M. A. Pfaller (ed.). Manual of clinical microbiology, 9<sup>th</sup>ed. American Society for Microbiology, Washington, D.C.
11. Mazura-Reetz, G. T. Neblett, and J. M. Galperin. 1979. MacConkey Agar: CO<sub>2</sub> vs. ambient incubation. Abst. Ann. Mtg. American Society for Microbiology. C179.
12. Isenberg, H. D. (ed.). 1992. Interpretation of aerobic bacterial growth on primary culture media, Clinical microbiology procedures handbook, vol.1, p. 1.6.1-1.6.7. American Society for Microbiology, Washington, D.C.

## AMBALAŽA/DOSTUPNOST

**BD MacConkey II Agar / BD Columbia CNA Agar Improved II with 5% Sheep Blood (Biplate)**

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© 2015 Becton, Dickinson and Company
**Becton Dickinson GmbH**