# **BD BACTEC™ Plus Aerobic/F Culture Vials** **Soybean-Casein Digest Broth in a Plastic Vial**

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## INTENDED USE

**BD BACTEC™ Plus Aerobic/F** medium is used in a qualitative procedure for the aerobic culture and recovery of microorganisms (bacteria and yeast) from blood. The principal use of this medium is with the **BD BACTEC** fluorescent series instruments.

## SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the **BD BACTEC** fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO<sub>2</sub> produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO<sub>2</sub> present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

Resins have been described for the treatment of blood specimens both prior to and after their inoculation into culture media. Resins have been incorporated into **BD BACTEC** culture media to enhance recovery of organisms without a need for special processing.<sup>1-3</sup>

## PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the **BD BACTEC** vial, CO<sub>2</sub> will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO<sub>2</sub> are monitored by the **BD BACTEC** fluorescent series instrument. Analysis of the rate and amount of CO<sub>2</sub> increase enables the **BD BACTEC** fluorescent series instrument to determine if the vial is positive, i.e., that the test sample contains viable organisms.

## REAGENTS

The **BD BACTEC** culture vials contain the following reactive ingredients prior to processing:

### List of Ingredients **BD BACTEC™ Plus Aerobic/F (442023)**

|                                   |              |
|-----------------------------------|--------------|
| Processed Water                   | 30 mL<br>w/v |
| Soybean-Casein Digest Broth       | 3.0%         |
| Yeast Extract                     | 0.25%        |
| Amino Acids                       | 0.05%        |
| Sugar                             | 0.2%         |
| Sodium Polyanetholsulfonate (SPS) | 0.05%        |
| Vitamins                          | 0.025%       |
| Antioxidants/Reductants           | 0.005%       |
| Nonionic Adsorbing Resin          | 13.4%        |
| Cationic Exchange Resin           | 0.9%         |

All **BD BACTEC** media are dispensed with added CO<sub>2</sub>.

## Warnings and Precautions:

The prepared culture vials are for *in vitro* diagnostic use.

This product contains dry natural rubber.

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens.**

**"Standard Precautions"<sup>4-7</sup> and institutional guidelines should be followed in handling all items contaminated with blood.**

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. DO NOT USE any vial showing evidence of contamination. A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, gas or contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. If a direct draw procedure is used, monitor the process closely to avoid refluxing materials into the patient.

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials that are cracked or leaking, or display turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions, a vial neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or **BD Luer-Lok™** tips.

#### Storage Instructions

The **BD BACTEC** vials are ready for use as received and require no reconstitution or dilution. Store in a cool, dry place (2–25 °C), **out of direct light**.

#### SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The recommended specimen volume is 8–10 mL. It is recommended that the specimen be inoculated into the **BD BACTEC** vials at bedside. A 10cc or 20cc syringe with a **BD Luer-Lok** brand tip is used to draw the sample, or a **BD Vacutainer®** Brand Needle Holder and a **BD Vacutainer** Brand Blood Collection Set, **BD Vacutainer Safety-Lok™** Blood Collection Set or other tubing “butterfly” set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. **The inoculated BD BACTEC vial should be transported to the laboratory as quickly as possible.**

#### PROCEDURE

Remove the flip-off cap from the **BD BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness in medium, and bulging or indented septums. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **NOT** recommended). Aseptically inject or draw directly 8–10 mL of specimen per vial. If sample volumes of 3–7 mL are used, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated aerobic vials should be placed in the BD BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the **BD BACTEC** fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive vial.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the **BD BACTEC** fluorescent series instrument and identified as such (see the appropriate **BD BACTEC** Fluorescent Series Instrument User’s Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the **BD BACTEC** fluorescent series instrument can determine a difference in fluorescence.

If at the end of the testing period a negative vial appears visually positive (i.e., chocolateized blood, bulging septum, lysed and/or very darkened blood in **BD BACTEC** Plus Aerobic/F medium), it should be subcultured and Gram-stained and treated as a presumptively positive vial.

Positive vials should be subcultured and Gram-stained. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to solid media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the **BD BACTEC** vials.

**Subculturing:** Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, use an appropriate venting unit (BD catalog # 249560) or equivalent. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

#### QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory’s standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

**DO NOT USE** culture vials past their expiration date.

**DO NOT USE** culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Standard M22, *Quality Control for Commercially Prepared Microbiological Culture Media*. The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

- |   |  |
|---|--|
| • <i>Neisseria meningitidis</i><br>ATCC 13090   | • <i>Candida glabrata</i><br>ATCC 66032      |
| • <i>Haemophilus influenzae</i><br>ATCC 19418   | • <i>Staphylococcus aureus</i><br>ATCC 25923 |
| • <i>Streptococcus pneumoniae*</i><br>ATCC 6305 | • <i>Escherichia coli</i><br>ATCC 25922      |
| • <i>Streptococcus pyogenes</i><br>ATCC 19615   | • <i>Alcaligenes faecalis</i><br>ATCC 8750   |
| • <i>Pseudomonas aeruginosa</i><br>ATCC 27853   |  |

\*CLSI recommended strain

For information on Quality Control for the **BD BACTEC** fluorescent series instrument, refer to the appropriate **BD BACTEC** Fluorescent Series Instrument User's Manual.

## RESULTS

A positive sample is determined by the **BD BACTEC** fluorescent series instrument and indicates the presumptive presence of viable microorganisms in the vial.

## LIMITATIONS OF THE PROCEDURE

### Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BD BACTEC** vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc.

### Recovery of SPS Sensitive Organisms From Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS, the presence of maximum volumes of blood (8–10 mL) can help to optimize recovery of these organisms. To enhance the growth of SPS sensitive organisms when less than 8 mL of blood is inoculated, additional whole human blood may be added.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is 3.0 mL or less, an appropriate supplement may be required for recovery of these organisms.

**BD BACTEC FOS™** Fastidious Organism Supplement or whole human blood may be used as nutritional supplements.

### Nonviable Organisms

A Gram-stained smear from a culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.

### Antimicrobial Activity

Neutralization of the antimicrobial activity by resins varies depending upon dosage level and timing of specimen collection. The use of supplementary additives should be considered in appropriate situations; as an example, the addition of penicillinase when  $\beta$ -lactam therapy is being employed.

### Recovery of *Streptococcus pneumoniae*

In aerobic media, *S. pneumoniae* will typically be visually and instrument positive, but in some cases no organism will be seen on Gram stain or recovered on routine subculture. If an anaerobic vial was also inoculated, the organism can usually be recovered by performing an aerobic subculture of the anaerobic vial, since this organism has been reported to grow well under anaerobic conditions.<sup>11</sup>

### General Considerations

Optimum recovery of organisms will be achieved by adding maximum amounts of blood. Published clinical studies have shown that the use of lower blood volumes may adversely affect recovery and/or detection times of organisms.<sup>12</sup> Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO<sub>2</sub> to be detected by the system or if significant growth has occurred before placing the vial into the system. In analytical studies, this device recovered 17 of 18 *Leuconostoc* spp. tested. False positivity may occur when the white blood cell count is high. The default 5 day protocol was utilized for all analytical testing with this device, and protocol lengths of >5 days have not been evaluated.

Due to the nature of biological materials in media products and inherent organism variability, the user should be cognizant of potential variable results in the recovery of certain microorganisms.

## EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

Performance of the **BD BACTEC** Plus Aerobic/F medium in glass vials has been established by a number of external clinical studies.<sup>1-3,8,9</sup> Seeded laboratory studies performed by BD have shown equivalent performance of the **BD BACTEC** Plus Aerobic/F medium in plastic vials compared to the **BD BACTEC** Plus Aerobic/F medium in glass vials.<sup>10</sup> The yeasts *Candida albicans*, *C. glabrata*, and *Cryptococcus neoformans* were tested in the analytical testing of this device.

## AVAILABILITY

### Cat. No. Description

442023 **BD BACTEC™** Plus Aerobic/F Medium

## REFERENCES

1. Wallis, C. et al. 1980. Rapid isolation of bacteria from septicemic patients by use of an antimicrobial agent removal device. *J. Clin. Microbiol.* 11:462-464.
2. Applebaum, P.C. et al. 1983. Enhanced detection of bacteremia with a new BACTEC resin blood culture medium. *J. Clin. Microbiol.* 17:48-51.
3. Pohlman, J.K. et al. 1995. Controlled clinical comparison of Isolator and BACTEC 9240 Aerobic/F resin bottle for detection of bloodstream infections. *J. Clin. Microbiol.* 33:2525-2529.
4. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.
6. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.

7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
8. Smith, J.A. et al. 1995. Comparison of BACTEC 9240 and BacT/Alert blood culture systems in an adult hospital. J. Clin. Microbiol. 33:1905-1908.
9. Flayhart, D. et al. 2007. Comparison of BACTEC Plus blood culture media to BacT/Alert FA blood culture media for detection of bacterial pathogens in samples containing therapeutic levels of antibiotics. J. Clin. Microbiol. 45:816-821.
10. Data available from BD Life Sciences.
11. Howden, R.J. 1976. Use of anaerobic culture for the improved isolation of *Streptococcus pneumoniae*. J. Clin. Pathol. 29:50-53.
12. Istrup, Duane M. and John A. Washington. 1983. The importance of volume of blood cultured in the detection of bacteremia and fungemia. Diagn. Microbiol. Infect Dis. 1:107-110.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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ATCC is a trademark of the American Type Culture Collection.
© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.
Benex Limited