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
1. Amies C.R. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, July 1967. Vol. 58, 296 - 300.
2. Isenberg H. D., Schoenkencht F.D. and Von Graeventiz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
3. Balows A., Hausler, Jr. W. J., Hermann K.L., Isenberg H. D., Shadomy H.J. Manual of Clinical Microbiology. Fifth Edition. American Society for Microbiology, Washington DC, 1991.
4. Isenberg H. D. (Editor in Chief). Clinical Microbiology Procedures Handbook. American Society for Microbiology, Washington DC, 1992.
5. Zavala M. K., Citron D. M., Goldstein E. J. C. Evaluation of a novel specimen transport system for anaerobic bacteria. Clinical Infectious Diseases, 1997; Vol. 25 (Suppl 2); S132 - 3.
6. Perry J.L. Assessment of swab transport systems for aerobic and anaerobic organism recovery. Journal of Clinical Microbiology, May 1997. Vol. 35, 1269 - 1271.
7. Van Horn K., Tóth C. and Wegienek J. Viability of aerobic microorganisms in four swab systems. Poster Session 249/C Abstract C-436. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
8. Van Horn K., Tóth C. and Warren K. Comparison of Copan Amies agar swab and BBL Port-a-cul swab for recovery of anaerobic bacteria. Poster Session 249/C Abstract C-437. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
9. Arbiq J.C, Forward K. R., and Le Blanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. Diagnostic Microbiology and Infectious Disease, 2000 Vol. 26, 163-168.
10. Van Horn K., Warren K., and Tóth C. Comparison of four swab systems for the recovery of anaerobic bacteria. Abstract C-32. 99th General Meeting of American Society for Microbiology, Chicago, IL, May 1999.
11. Olsen C. C., Schwabke J. R., Benjamin W. H., Beverly A., and Waites K. B. Comparison of direct inoculation and Copan transport systems for the isolation of *Neisseria gonorrhoeae* from endocervical specimens. Journal of Clinical Microbiology Vol. 37, No. 11, Nov. 1999, p. 3583 - 3585.
12. Citron D.M., Warren Y.A., Hudspeth M.K. and Goldstein E.J.C. Survival of aerobic and anaerobic bacteria in purulent clinical specimens maintained in Copan Venturi Transystem and Becton Dickinson Port-a-Cul transport systems. Journal of Clinical Microbiology Vol. 38, No. 2. Feb. 2000, p.892-894.
13. Thompson D. S., French S. A. Comparison of commercial Amies transport systems with in-house Amies medium for recovery of *Neisseria gonorrhoeae*. Journal of Clinical Microbiology Vol. 37, No. 9. Sept. 1999, p. 3020-3021.
14. Hetchler C., Brown C., and Galbraith J. Comparison of 3 Amies Gel transport systems for the recovery of 12 clinically significant organisms. Abstract C-152. 100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.
15. Hindiyeh M., Acevedo V., Croft A., and Carroll K. Comparison of the new Starplex StarSwab II and the new Copan Vi-Pak Amies Agar Gel collection and transport swabs with BBL Port-a-Cul for the maintenance of anaerobic and fastidious aerobic organisms. Abstract C-154. 100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.
16. Summanen P., Baron E.J., Citron D. M., Strong C.A., Wexler H.M. and Finegold S. M. Wadsworth Anaerobic Bacteriology Manual. Fifth Edition. Star Publishing Company, Belmont, CA, 1993.

© 2017 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.

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**R<sub>x</sub> Only** This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

**CE 0123** Identification number of notified body / Identification de l'organisme notifié / Identifizierung der benannten Stelle / Identificación del organismo notificado / Identificazione dell'organismo notificato / Identificação do organismo notificado / Identifikační číslo notifikovaného subjektu / Identifikasjonsnummer ved påvist organisme / Identifieringsnummer av certifierad myndighet

**STERILE R** Method of sterilization: irradiation / Метод на стерилизация: ирадиация / Způsob sterilizace: záření / Steriliseringmetode: bestråling / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Método de esterilización: irradiación / Steriliseerimismeetod: kiirgus / Méthode de stérilisation : irradiation / Metoda sterilizacije: zračenje / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Стерилизация әдісі – сәуле түсіру / 소독 방법: 방사 / Sterilizavimo būdas: radiacija / Sterilizēšanas metode: apstarošana / Gesteriliseerd met behulp van bestraling / Steriliseringmetode: bestråling / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metodă de sterilizare: iradiere / Метод стерилизации: облучение / Metóda sterilizácie: ožiarenie / Metoda sterilizacije: ozračavanje / Steriliseringmetod: strålning / Sterilizasyon yöntemi: iradyasyon / Метод стерилизації: опромінення / 灭菌方法: 辐射

**⌘** Do not reuse / Не използвайте отново / Nepoužívejte opakovaně / Ikke til genbrug / Nicht wiederverwenden / Μην επαναχρησιμοποιείτε / No reutilizar / Mitte kasutada korduvalt / Ne pas réutiliser / Ne koristiti ponovo / Egyszer használatos / Non riutilizzare / Пайдаланбаңыз / 재사용 금지 / Tik vienkartiniam naudojimui / Nelietot atkārtoti / Niet opnieuw gebruiken / Kun til engangsbruk / Nie stosować powtórnie / Não reutilize / Nu refolositi / Не использовать повторно / Nepoužívejte opakovane / Ne upotrebľavajte ponovo / Får ej återanvändas / Tekrar kullannayin / Не використовувати повторно / 请勿重复使用

**🌡** Temperature limitation / Температурни ограничения / Teplotní omezení / Temperaturbegrænsning / Temperaturbegrenzung / Περιορισμοί θερμοκρασίας / Limitación de temperatura / Temperatuuri piirang / Limites de température / Dozvoljena temperatura / Hőmérsékleti határ / Limiti di temperatura / Температураны шектеу / 온도 제한 / Laikymo temperatūra / Temperatūras ierobežojumi / Temperatuurlimit / Temperaturbegrensning / Ograniczenie temperatury / Limites de temperatura / Limite de temperatură / Ограничение температуры / Ohraničenje teploty / Ograničenje temperature / Temperaturgräns / Sıcaklık sınırlaması / Обмеження температури / 温度限制

**👉** Peel / Обелете / Otevfete zde / Åbn / Abziehen / Αποκολλήστε / Desprender / Koorida / Décoller / Otvoriti skini / Húzza le / Staccare / Ўстиґи қабатын алып таста / 벗기기 / Plěšti čia / Atliimēt / Schillen / Trekk av / Oderwać / Destacar / Se dezlipeste / Отклеить / Odtrhnite / Oljuštiti / Dra isār / Ayırma / Відклеїти / 撕下

**Manufactured by:**  
Copan Italia SpA  
Via F. Perotti, 10  
25125 Brescia Italy

**Distributed by:**  
Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA

**Australian Sponsor:**  
Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia

ATCC is a trademark of the American Type Culture Collection.

© 2017 BD, BD, and the BD Logo are trademarks of Becton, Dickinson and Company.