

**PID**

**NAMJENA**

**BD Phoenix** PID panel namijenjen je za brzo otkrivanje većine aerobnih i fakultativnih anaerobnih gram pozitivnih bakterija humanog porijekla (za potpuni popis uzoraka pogledajte Korisnički priručnik Sustava **BD Phoenix**). Ovaj panel namijenjen je za uporabu samo s instrumentom **BD Phoenix** automatizirani mikrobiološki sustav.

**NAČELA POSTUPKA**

Pogledajte Korisnički priručnik Sustava **BD Phoenix**.

**PANELI**

Pogledajte Korisnički priručnik Sustava **BD Phoenix**.

**POHRANA**

Pohraniti na 15–25 °C. Ne koristiti panel ako je vrećica probušena ili otvorena. Ne rabiti panel ako nedostaje sredstvo za sušenje ili ako je vrećica sa sredstvom za sušenje potrgana.

**PANELI SE MORAJU ISKORISTITI U ROKU OD 2 H NAKON VAĐENJA IZ VREĆICE.**

**Mjere opreza:** Za *in vitro* dijagnostiku.

**UPOZORENJA**

U svim postupcima poštujujte uobičajene mjere opreza za mikrobiološke opasnosti. Pri rukovanju svim predmetima kontaminiranim uzorcima i mikroorganizmima treba se pridržavati „Standardnih mjera opreza“<sup>1,2</sup> i institucionalnih smjernica. Prije bacanja sterilizirajte spremnike uzoraka te druge inokulirane materijale autoklavom.

**KONTROLA KVALITETE**

U nastavku pogledajte očekivane rezultate.

Ispitni organizam	Očekivani rezultat
<i>Staphylococcus aureus</i> ATCC 29213	<i>Staphylococcus aureus</i>
<i>Enterococcus faecalis</i> ATCC 29212	<i>Enterococcus faecalis</i>

**FUNKCIONALNE KARAKTERISTIKE**

Pogledajte Korisnički priručnik Sustava **BD Phoenix**.

**REFERENCE**

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. U.S. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. 1991. Occupational exposure to bloodborne pathogens. Federal Register 56:64175–64182.

Tehnički servis i podrška: obratite se lokalnom predstavniku tvrtke BD ili posjetite [www.bd.com](http://www.bd.com).

ATCC is a trademark of the American Type Culture Collection.

© 2016 BD. BD, the BD Logo and BD Phoenix are trademarks of Becton, Dickinson and Company.



Becton, Dickinson and Company  
7 Loveton Circle  
Sparks, MD 21152 USA



Benex Limited  
Pottery Road, Dun Laoghaire  
Co. Dublin, Ireland

**Australian Sponsor:**

Becton Dickinson Pty Ltd.  
4 Research Park Drive  
Macquarie University Research Park  
North Ryde, NSW 2113  
Australia



MADE IN USA  
L006814(04)  
2016-05